Rainbow™ Mini-Loop-Fitnessbänder

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Rainbow™ Latex-Free Mini Loop

Intended Purposes Wrap around upper or lower extremities for rehabilitation, activation, and strengthening of muscles.

Indications Use for physical rehabilitation, muscle strengthening and activation, or range of motion improvement.

Contraindications Resistance training may not be recommended for those with cardiac issues or hypertension.

Instructions for Use Place Mini Loop exercise band around body part. Flatten the band for the most comfort. Flex or extend body part against resistance, while maintaining balance and control.

Instructions for Care Hand wash with mild soap, dry flat.

Warning

- Exercise bands and tubing should be used for their intended use only.
 - Before beginning this or any type of exercise program, please consult your health care professional.
 - Not to be used by unsupervised children.
 - Before each exercise session, inspect the resistance band for small nicks, tears, or punctures that may cause the band to break. If you find any flaws, discard the product and replace before performing any exercises.
 - Do not wear jewelry that may tear or damage band.
 - Use caution with long or sharp fingernails.
 - Do not wrap band/tubing around uneven surfaces.
 - Do not stretch exercise band/tubing more than two times normal length. Doing so may cause the band/tubing to snap.
 - Snapping bands/tubing can cause serious injury, including fractures, eye, nerve and muscle damage.
 - Wearing protective eyewear is recommended.
 - As an added precaution, keep band/tubing and accessories away from head and face while exercising.
 - When using with an accessory anchored to a door, make sure that the door opens away from you.
 - Only exercise on smooth and debris free surfaces.
 - Ensure that bands are securely fitted and/or anchored to accessories.
 - As with any exercise program, you may experience muscle soreness over the first few days. If your pain persists, stop your exercise program and consult with health care professional.
 - Do not exercise while experiencing pain.
 - Store your exercise band/tubing in a protective bag or box.
 - Do not expose band/tubing to direct heat or sunlight.
 - Detach band/tubing from exercise handles and accessories before storing.
- With heavy use, such as in a physical therapy clinic, band/tubing should be replaced every 1-2 months.

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North Coast Medical, Inc.
Morgan Hill, CA 95037 – U.S.A.
ICON (LR) Limited
Allenspach Medical AG