



North Coast

Comfort Gel Pads



Thermoplastic splint with adhesive-backed Curved Pad



Comfort Cool™ splint with hook-backed Rectangular Pad

Instructions For Use

1. Do not apply the Comfort Gel Pad over an open wound. The pad can be used when the wound is closed, dry and sutures have been removed.
2. Wash hands and the involved area in mild, soapy water. Pat dry.
3. Apply **adhesive back pad** to molded thermoplastic splint in alignment with the scar or vulnerable area. Do not immerse in hot water as this can cause the adhesive backing to release. Once placed, the adhesive is difficult to relocate. (Figure 1)
4. Apply **hook-back pad** to the fabric lining of prefabricated splints in alignment with the scar or vulnerable area. The pad can be repositioned as needed (Figure 2).
5. If needed, the gel pad can be trimmed with scissors.

Instructions For Care

Washing the surface of the pad will reactivate the tackiness for better adherence to the skin if the pad becomes overly embedded with surface dirt.

Hand wash the skin area and the pad daily for best results. Use mild soap in cool water and air dry.

Rectangular Pad

NC85050A	Adhesive back
NC85050H	Hook-back

Curved Pad

NC85051A	Adhesive back
NC85051H	Hook-back

To be used under the guidance of a qualified medical professional.



North Coast Medical, Inc.
Gilroy, California U.S.A.
Authorised Representative
Medica Surgical Innovations Ltd
BB2 4PB UK