



Comfortprene™ Gel Support


Indications:

Use to provide moderate support and compression around the palm, wrist and thumb CMC. Separate gel pad provides protective padding and scar care anywhere. Ideal for post Carpal Tunnel Release, de Quervain's Release and after other surgical procedures along the radial thumb or hand. The gel pad is made of hypoallergenic silicone gel with hook backing. The support is made of 1/8" (3.2mm) neoprene with terry liner. Contains no latex.



Instructions For Use:

1. Do **not** apply the gel pad over an open wound. The pad can be used when the wound is closed, dry and sutures have been removed. Wash hands and the involved area in mild, soapy water.



2. Apply hook-back pad to the fabric lining of the support in alignment with the scar or vulnerable area. (Figure 1). The pad can be repositioned and trimmed with scissors if needed. Once placed correctly, the gel pad can be used up to six weeks.
3. To apply the Comfortprene™ support, place the wider mid section over the palm and the thinner extension through the web space to the back of the hand. Bring the radial section around to the back of the hand and fasten it to the hook at the end of the web space section. Bring the ulnar section around and fasten it to the hook strap on the back of the hand.

Care Instructions:

1. Wash the skin and gel pad daily for best results. Use mild soap, cool water and air dry. Washing the pad's surface will reactivate tackiness for better adherence to the skin until the gel pad becomes too embedded with surface dirt.
2. Hand wash the support in cool, soapy water and air dry.

Comfortprene™ Gel Support: Beige

Left	Right	
NC22641-1	NC22641-2	Small
NC22641-3	NC22641-4	Medium
NC22641-5	NC22641-6	Large

Comfortprene™ Gel Support: Black

Left	Right	
NC22642-1	NC22642-2	Small
NC22642-3	NC22642-4	Medium
NC22642-5	NC22642-6	Large

Comfortprene™ Gel Support: Blue

Left	Right	
NC22643-1	NC22643-2	Small
NC22643-3	NC22643-4	Medium
NC22643-5	NC22643-6	Large

NC85050H Hook Replacement Gel Pad

To be used under the guidance of a qualified medical professional.



North Coast Medical, Inc.
Gilroy, California U.S.A.
Authorised Representative
Medica Surgical Innovations Ltd.
BB2 4PB UK