





# Line Stoppers

## Indications

Use to attach monofilament lines to slings, loops, Gyovai finger springs and rubber bands. Allows the line to be locked into place without using knots. The line cannot slide when tension is applied in one direction but length can still be adjusted. Made of aircraft quality aluminum.

## Instructions For Assembly

1. Determine the number of lines to use based on the information listed below.  
In general, finger slings have 4 lines, finger loops have 2 lines, and Gyovai finger spring's line is looped back to make 2 lines.
2. With the lines together, weave the monofilament line(s) up and through the 2 smaller holes in the metal piece.
3. Check the locking hold by pulling on the line(s) leading off the small hole end of the piece, opposite the big hole end. These lines are usually attached to the finger sling/loop, or springs.
4. The unattached end of the line extending out from the Line Stopper can be cut down. However, this is the line to use if the distance to the sling/loop needs to be lengthened.
5. To adjust the line length, simply push the line(s) toward the Line Stopper and lift the line(s) between the 2 small holes, then pull the end of the monofilament line(s) leading to or away from the attachment sling/loop.
6. Use the Line Stopper's big hole to attach rubber bands or elastic cord via regular, square knot or slipknot securement.



**Line Recommendations** — Holes are sized to manage certain size monofilament lines for better locking results.

**NC22556 Small Hole Line Stopper** – Use with 10-40 lb monofilament line.

Monofilament Size	Number of lines through the holes
10 lb	4x
15 lb	2x
20 lb	2x
30 lb	1x
40 lb	1x

<b>Suggested Monofilament Lines</b>	
<i>Nylon Outrigger Line</i>	
NC12511	10lbs line
NC12511-15	15lbs line
NC12511-20	20lbs line
<i>Gyovai Finger Spring</i>	
NC12508	15lbs line

**NC22557 Large Hole Line Stopper** – Use with 20-40 lb monofilament line.

Monofilament Size	Number of lines through the holes.
20 lb	2x
30 lb	1x
30 lb	2x
40 lb	1x
40 lb	2x

To be used under the guidance of a qualified medical professional.



North Coast Medical, Inc.  
Gilroy, California U.S.A.  
Authorised Representative  
Medica Surgical Innovations Ltd.  
BB2 4PB UK