



North Coast

Brands that Perform™

### Freehand Wrist Brace®

*Adjustable, low-profile dorsal splint for treating carpal tunnel syndrome.*

The Freehand Wrist Brace® is designed to prevent extreme wrist flexion and help maintain a neutral wrist, while allowing full hand use without contact or pressure over the carpal tunnel. The soft tubing gently fits around the middle finger and attaches to a 1" (2.5cm) wide dorsal leather-bound forearm support. The soft hook and loop straps can be positioned anywhere on the support piece for proper fit and comfort. The open design of the wrist brace allows the hand to reach into pockets or other tight places while also keeping the brace cool and dry. The Freehand Wrist Brace® comes in one universal size which can fit most hands and is hand washable without the need to remove the splint. Latex-free.

HW465 15.40



### Freehand CMZ Brace®

*Provides an adjustable, unobtrusive fit to support the thumb CMC joint.*

The Freehand CMZ Brace® positions the thumb CMC joint into abduction for such conditions as osteoarthritis or ligament injuries. The preformed, bendable wire encased in tubing provides maximum comfort due to its flexibility, springiness, and lack of pressure over the CMC joint. With its ease of application, removal, and adjustment the brace can be easily modified for a customized fit. The Freehand CMZ Brace® comes in one universal size which can fit most hands and is hand washable without the need to remove the splint. Latex-free.

HW467 15.40



### Freehand Dex Dynamic Elbow Extension Brace

*Low-profile, lightweight, dynamic support can be used for a variety of elbow conditions.*

The Freehand Dex Dynamic Elbow Extension Brace is ideal for treating those with cubital tunnel syndrome. The brace prevents prolonged elbow flexion and is specifically designed to avoid any pressure on the ulnar nerve. The side coils provide dynamic extension force allowing the temporarily flexed elbow to return to extension, even during sleep. This brace allows momentary elbow flexion against slight resistance for functional arm use, which can also help maintain or gently strengthen the biceps.

The coils and side pieces are made of preformed bendable wire encased in tubing. The leather covered moldable crossbars allow for an easy, custom fit that gently contours around the arm. The open design results in minimal skin contact allowing the arm to remain cool and dry. The Freehand Dex Dynamic Elbow Extension Brace comes in two sizes to fit arms of all lengths and girths. Fits either the left or right elbow. Latex-free.

HW469-S/M 37.95

HW469 L/XL 37.95