

**New**

*OrthoThermic™ Thumb CMC  
Restriction Splint*



**Ortho  
Thermic™**



***Soothing  
Comfort***

**COMPRESSION + HEAT**

# Ortho Thermic™

THERAPEUTIC  
PAIN RELIEF

## Thumb CMC Restriction Splint



## Soothing Comfort

Provides support and warmth for the thumb CMC joint.

Features a moldable metal stay for enhanced support  
at the thumb's base.

Silicon beads on outer glove surface improve gripping ability.  
Heat-generating antimicrobial liner.

Not made  
with natural  
rubber latex.

OrthoThermic™ Thumb CMC Restriction Splint

Left	Right	Size	MP Circumference	
NC79670	NC79671	Small	6" - 7"	(15-18 cm)
NC79672	NC79673	Small +	6½" - 7¼"	(17-18 cm)
NC79674	NC79675	Medium	7" - 8"	(18-20 cm)
NC79676	NC79677	Medium +	7¼" - 8¼"	(18-21 cm)
NC79678	NC79679	Large	8" - 9"	(20-23 cm)
NC79680	NC79681	Large +	8¾" - 9½"	(21-24 cm)
NC79682	NC79683	X-Large	9" - 10"	(23-25 cm)



To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call 800.821.9319

Available in Virtual Showroom™

COMPRESSION + HEAT

# OrthoThermic™

THERAPEUTIC  
PAIN RELIEF

## Gloves Soothing Comfort

COMPRESSION + HEAT



Provides a thermal effect to increase body heat.

Wicks away moisture.

Silicon beads on outer glove surface for improved gripping ability.

Antimicrobial liner.

Sold in pairs.

Hand wash with mild soap and water. Air dry.



SCAN  
to learn more



OrthoThermic™ Gloves  
Size

MP Circumference

	Size	MP Circumference	
NC54100-XS	X-Small	6" - 6¾"	(15-17 cm)
NC54100-SM	Small	7" - 7¾"	(18-20 cm)
NC54100-MD	Medium	8" - 9"	(21-23 cm)
NC54100-LG	Large	9¼" - 10½"	(24-26 cm)
NC54100-XL	X-Large	10¾" - 11½"	(27-29 cm)
NC54100-XXL	XX-Large	Over 11¾"	(Over 30 cm)

To order, log in at [www.ncmedical.com](http://www.ncmedical.com)  
or call **800.821.9319**

Available in Virtual Showroom™

# Compression Gloves

**Provides light to moderate compression for the hand and wrist.**

- Provides warmth and soft, gentle compression to help soothe soft tissue edema and arthritic joint pain.
- Two wrist lengths and two finger styles feature an open slit in the cuff, for easy application and non-constricting fit.
- Made of a soft, silky nylon and spandex blend, with seams on the outside for smooth edema control.
- Sold individually, Left or Right.
- Machine wash delicate, air dry.

**Now available in *XL* sizes**



**Beige  
Full Finger  
Over The Wrist Length**

	Left	Right
XS	NC53210	NC53211
S	NC53212	NC53213
M	NC53214	NC53215
L	NC53216	NC53217
XL	NC53218	NC53219

**Beige  
Full Finger Wrist Length**

	Left	Right
XS	NC53206	NC53207
S	NC53200	NC53201
M	NC53202	NC53203
L	NC53204	NC53205
XL	NC53208	NC53209

**Sizing chart  
available online.**

Not made with  
natural rubber latex.



**Black  
Tipless Over The Wrist Length**

	Left	Right
XS	NC53220-B	NC53221-B
S	NC53222-B	NC53223-B
M	NC53224-B	NC53225-B
L	NC53226-B	NC53227-B
XL	NC53228-B	NC53229-B

**Beige  
Tipless Over The Wrist Length**

	Left	Right
XS	NC53220	NC53221
S	NC53222	NC53223
M	NC53224	NC53225
L	NC53226	NC53227
XL	NC53228	NC53229

**To order, log in at [www.ncmedical.com](http://www.ncmedical.com)  
or call **800.821.9319****

Available in Virtual Showroom™

# Heavy Compression Gloves

**Moderate to heavy compression  
for the hand and wrist.**

- Outside seams prevent unwanted pressure lines.
- Open slit in the cuff for easy application and non-constricting fit.
- Made of 83% nylon and 17% spandex.
- Wrist length measures 3" (7.6 cm) proximal to wrist crease.

*Available individually and in pairs.*

## Full Finger - Singles

	Size	Hand	MP Circumference
NC55320	X-Small	Left	Less than 7" (Less than 18 cm)
NC55321	X-Small	Right	Less than 7" (Less than 18 cm)
NC55322	Small	Left	7" to 8" (18 to 20 cm)
NC55323	Small	Right	7" to 8" (18 to 20 cm)
NC55324	Medium	Left	7½" to 8½" (19 to 22 cm)
NC55325	Medium	Right	7½" to 8½" (19 to 22 cm)
NC55326	Large	Left	8" to 10" (20 to 25 cm)
NC55327	Large	Right	8" to 10" (20 to 25 cm)

## Open Finger - Singles

	Size	Hand	MP Circumference
NC55310	X-Small	Left	Less than 7" (Less than 18 cm)
NC55311	X-Small	Right	Less than 7" (Less than 18 cm)
NC55312	Small	Left	7" to 8" (18 to 20 cm)
NC55313	Small	Right	7" to 8" (18 to 20 cm)
NC55314	Medium	Left	7½" to 8½" (19 to 22 cm)
NC55315	Medium	Right	7½" to 8½" (19 to 22 cm)
NC55316	Large	Left	8" to 10" (20 to 25 cm)
NC55317	Large	Right	8" to 10" (20 to 25 cm)

## Full Finger - Pairs

	Size
NC55331-XS	X-Small
NC55331-SM	Small
NC55331-MD	Medium
NC55331-LG	Large

## Open Finger - Pairs

	Size
NC55330-XS	X-Small
NC55330-SM	Small
NC55330-MD	Medium
NC55330-LG	Large

*Machine washable.*

*Not made with  
natural rubber latex.*



To order, log in at **www.ncmedical.com**  
or call **800.821.9319**

Available in Virtual Showroom™



# Thumb CMC Restriction Splint

Lightweight, comfortable support.

All Comfort Cool® products are made of a unique, perforated neoprene material lined with soft terry cloth featuring anti-microbial technology.

Comfort Cool® products are an innovative combination of design and materials that make them the best choice for comfort and compliance.

The perforated neoprene material offers soft cushioning for comfort; perforations provide ventilation for more breathability while the terry cloth lining helps keep skin cool.

## BLACK

Also available in  
Pediatric sizes.



Comfort Cool® Thumb CMC Restriction, BLACK

	Size	DPC Circumference	Hand
NC79562	Sm	6" to 7" (15 to 18 cm)	Left
NC79563	Sm	6" to 7" (15 to 18 cm)	Right
NC79556	Sm Plus	6½" to 7¼" (17 to 18 cm)	Left
NC79557	Sm Plus	6½" to 7¼" (17 to 18 cm)	Right
NC79564	Med	7" to 8" (18 to 20 cm)	Left
NC79565	Med	7" to 8" (18 to 20 cm)	Right
NC79558	Med Plus	7¼" to 8¼" (19 to 21 cm)	Left
NC79559	Med Plus	7¼" to 8¼" (19 to 21 cm)	Right
NC79566	Lg	8" to 9" (20 to 23 cm)	Left
NC79567	Lg	8" to 9" (20 to 23 cm)	Right
NC79560	Lg Plus	8¼" to 9½" (21 to 24 cm)	Left
NC79561	Lg Plus	8¼" to 9½" (21 to 24 cm)	Right
NC79550	X-Lg	9" to 10" (23 to 25 cm)	Left
NC79551	X-Lg	9" to 10" (23 to 25 cm)	Right

Provides direct support for the thumb CMC joint while allowing full finger function.

- Patented contoured neoprene strap wraps around the thumb CMC joint and the first metacarpal for direct CMC joint support.
- Specifically designed for osteoarthritis and medial subluxation of the thumb CMC joint.
- Reinforced stitching allows the splint to be trimmed below the thumb MP joint if desired.
- Provides gentle compression that helps relieve thumb discomfort associated with CMC joint arthritis and/or instability, tendinitis or repetitive motion.
- Made of thin, 1/16" (1.6 mm) perforated neoprene with a terry cloth liner that helps keep skin cool and comfortable.
- Plus sizes are a half size bigger in the hand and a full size larger in thumb circumference to accommodate enlarged thumb conditions.
- Not made with natural rubber latex.



Comfort Cool® Thumb CMC Restriction, BEIGE

	Size	DPC Circumference	Hand
NC79580	Sm	6" to 7" (15 to 18 cm)	Left
NC79581	Sm	6" to 7" (15 to 18 cm)	Right
NC79582	Sm Plus	6½" to 7¼" (17 to 18 cm)	Left
NC79583	Sm Plus	6½" to 7¼" (17 to 18 cm)	Right
NC79584	Med	7" to 8" (18 to 20 cm)	Left
NC79585	Med	7" to 8" (18 to 20 cm)	Right
NC79586	Med Plus	7¼" to 8¼" (19 to 21 cm)	Left
NC79587	Med Plus	7¼" to 8¼" (19 to 21 cm)	Right
NC79588	Lg	8" to 9" (20 to 23 cm)	Left
NC79589	Lg	8" to 9" (20 to 23 cm)	Right
NC79590	Lg Plus	8¼" to 9½" (21 to 24 cm)	Left
NC79591	Lg Plus	8¼" to 9½" (21 to 24 cm)	Right
NC79592	X-Lg	9" to 10" (23 to 25 cm)	Left
NC79593	X-Lg	9" to 10" (23 to 25 cm)	Right

US PATENT# 6,702,772

To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call 800.821.9319

Available in Virtual Showroom™



Lightweight, comfortable support.

# Wrist & Thumb CMC Restriction Splint

Wide neoprene strap provides extra support and compression around the wrist.



**Special pull-through wrist strap**  
*Requires less hand manipulation,  
making it easy to adjust and  
fasten the splint.*

- Supports the wrist for conditions such as tendinitis, arthritis, de Quervain's, joint instabilities, and post-surgical repair of the wrist and thumb.
- Reinforced stitching allows the splint to be trimmed below the thumb MP joint, if desired.
- Made of thin, 1/16" (1.6 mm) perforated neoprene with a terry cloth liner that helps keep skin cool and comfortable.
- Not made with natural rubber latex.

*Direct support for the thumb CMC joint  
while allowing full finger function.*



**Comfort Cool® Wrist and Thumb CMC Restriction Orthosis**

	Size	MP Circumference	Hand
NC79570	Sm	6" to 7" (15 to 18 cm)	Left
NC79571	Sm	6" to 7" (15 to 18 cm)	Right
NC79572	Med	7" to 8" (18 to 20 cm)	Left
NC79573	Med	7" to 8" (18 to 20 cm)	Right
NC79574	Lg	8" to 9" (20 to 23 cm)	Left
NC79575	Lg	8" to 9" (20 to 23 cm)	Right
NC79576	X-Lg	9" to 10" (23 to 25 cm)	Left
NC79577	X-Lg	9" to 10" (23 to 25 cm)	Right

To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call **800.821.9319**

Available in Virtual Showroom™



# Thumb CMC Abduction Orthosis

Lightweight, comfortable support.

All Comfort Cool® products are made of a unique, perforated neoprene material lined with soft terry cloth featuring anti-microbial technology.

Comfort Cool® products are an innovative combination of design and materials that make them the best choice for comfort and compliance.

The perforated neoprene material offers soft cushioning for comfort; perforations provide ventilation for more breathability while the terry cloth lining helps keep skin cool.

**Provides direct support for the thumb CMC joint while allowing full finger function.**

- Same orthosis as the Comfort Cool® Thumb CMC Restriction Splint, but with the addition of a firm, moldable insert inside the thumb strap.
- Ideal for treating thumb deformity and contractures due to neurological conditions, median nerve injuries or arthritis, especially where the thumb is held abducted against the palm.
- Contoured neoprene strap helps to abduct and stabilize the first metacarpal and open the thumb web space.
- Reinforced stitching allows the splint to be trimmed below the thumb MP joint if desired.
- Made of thin, 1/16" (1.6 mm) perforated neoprene with a terry cloth liner that helps keep skin cool and comfortable.
- Plus sizes are a half size bigger in the hand and a full size larger in thumb circumference to accommodate enlarged thumb conditions.
- Not made with natural rubber latex.



**Moldable insert**  
*Shape to fit precisely along the thenar eminence crease.*

**Comfort Cool® Thumb CMC Abduction Orthosis**

	Size	MP Circumference	Hand
NC79494	Sm	6" to 7" (15 to 18 cm)	Left
NC79495	Sm	6" to 7" (15 to 18 cm)	Right
NC79494-P	Sm Plus	6½" to 7¼" (17 to 18 cm)	Left
NC79495-P	Sm Plus	6½" to 7¼" (17 to 18 cm)	Right
NC79496	Med	7" to 8" (18 to 20 cm)	Left
NC79497	Med	7" to 8" (18 to 20 cm)	Right

	Size	MP Circumference	Hand
NC79496-P	Med Plus	7¼" to 8¼" (19 to 21 cm)	Left
NC79497-P	Med Plus	7¼" to 8¼" (19 to 21 cm)	Right
NC79498	Lg	8" to 9" (20 to 23 cm)	Left
NC79499	Lg	8" to 9" (20 to 23 cm)	Right
NC79498-P	Lg Plus	8¼" to 9½" (21 to 24 cm)	Left
NC79499-P	Lg Plus	8¼" to 9½" (21 to 24 cm)	Right

To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call 800.821.9319

Available in Virtual Showroom™

# Nite-Nite™ Neutral Wrist Support

**Provides a comfortable, neutral resting position for the wrist and hand.**

*Adjustable, removable metal stay helps to immobilize the wrist*



- Breathable, nylon-lined palm pillow stays cool for comfort while sleeping.
- Padding below wrist and forearm for added comfort.
- Soft polyester lining extends from forearm to the palm.
- Three closures for an adjustable fit.
- Ventilated to assist in releasing heat and moisture.
- Fits both left or right hand.
- One size fits most.
- Not made with natural rubber latex.

NC12220 Nite-Nite™ Neutral Wrist Support

*Breathable palm pillow*



**ALL NIGHT  
COMFORT AND  
PROTECTION**

To order, log in at **[www.ncmedical.com](http://www.ncmedical.com)** or call **800.821.9319**

Available in Virtual Showroom™

# Finger Extension Tube

**Reduce flexion contractures.**

- The shape and design of the Norco® Finger Extension Tube provides a low stress, high duration capability, ideal for flexion contractures.
- Circumferential pressure helps restrict fluid to aid in edema reduction.
- Tip of the finger should be visible. If needed, trim either end.
- Orthotic should fit comfortably, a firm fit is not required for effective results.
- Not recommended for flexion contractures greater than 50°.
- Made of dense neoprene, polyester, and poly/cotton tape.
- Not made with natural rubber latex.



*Poly/cotton tape  
provides comfort and durability.*

	Size	PIP Joint Circumference	
NC79068-1	X-Small	1¾" to 2"	(4.5 to 5.1 cm)
NC79068-2	Small	2" to 2½"	(5.1 to 6.4 cm)
NC79068-3	Medium	2½" to 2¾"	(6.4 to 7.0 cm)
NC79068-4	Red	12¾" to 3"	(7.0 to 7.6 cm)
NC79068-5	Red	3" to 3½"	(7.6 to 8.9 cm)

To order, log in at **[www.ncmedical.com](http://www.ncmedical.com)** or call **800.821.9319**

Available in Virtual Showroom™

# HOOK SPLINT.®

A sleek, low profile Relative Motion Orthosis.  
Relative motion made extremely simple.



Secure a custom fit with  
adjustable locking tabs.



Relative motion is a technique which places one or more fingers in a relatively extended or flexed position compared to the others.

Allows the entire hand to stay functional while achieving treatment objectives.

Formed by a metallic band with a smooth silicone finish, and a locking tab on each end.

Designed as a ready-to-use alternative to thermoplastics.

Waterproof, UV coating protects the outer surface and ensures long-term patient comfort.

#### Indications include:

Extensor tendon/sagittal band repair, PIP extension or flexion contracture, boutonniere deformity, trigger finger, Dupuytren's disease, tendon adhesions, digit pain/stiffness.



SCAN  
for VIDEO

	Size	Length	Width
<b>Singles</b>			
NC58774-1	SM/MED	11" (28 cm)	0.4" (1 cm)
NC58775-1	LG	11.8" (30 cm)	0.5" (1.2 cm)
NC58776-1	X-LG	14" (36 cm)	0.5" (1.2 cm)
NC58777-1	MAX	17.3" (44 cm)	0.5" (1.2 cm)
<b>Case of 10</b>			
NC58774-10	SM/MED	11" (28 cm)	0.4" (1 cm)
NC58775-10	LG	11.8" (30 cm)	0.5" (1.2 cm)
NC58776-10	X-LG	14" (36 cm)	0.5" (1.2 cm)
NC58777-10	MAX	17.3" (44 cm)	0.5" (1.2 cm)



To order, log in at **www.ncmedical.com**  
or call **800.821.9319**

Available in Virtual Showroom™

# ULTRA- ORTHO™

**Multipurpose Orthotic**



Provides wrist positioning and adjustability for both flexion and extension of the hand, wrist, thumb and digits.

Allows Metacarpophalangeal (MCP) joints to be placed into **extension** for extensor tendon injuries, or **flexion** for joint injuries.



C-BAR  
Maintains webspace of thumb allowing a good sweeping position for adjustable abduction, extension and flexion.



*"The Progress Ultra-Ortho is a functional solution for therapists and physicians. The orthotic is versatile and customizable for a wide variety of diagnoses, both orthopedic and neurological. The skilled craftsmanship and superior materials provide confidence knowing patients are getting a product that will therapeutically protect, progress, and meet all functional goals."*

- Greg Pitts, OTD, MSOTR/L, CHT



The innovative design combines a malleable frame with dense foam padding for a lightweight, soft orthosis that supports like a moldable plastic splint.

**ULTRA-  
ORTHO™**  
Multipurpose Orthotic



Independent thumb strap allows custom adjustability for securing and supporting thumb in place.

All straps are fastened with a rivet to meet L-Code requirements. Each strap is rounded to prevent shearing and offer maximum comfort and control.



Ultra-Ortho™ is biomechanically designed to immobilize and offer adjustable support to meet a variety of therapeutic goals.

The orthotic is easily adjustable and does not cause any impingement or pressure areas when straps are fastened. This allows the orthotic to be used immediately for a variety of different types of pathologies, maximizing the functional outcome and progress as needed without refitting for a new orthotic.

Applications include but are not limited to: contractures due to Dupuytren, head injury, stroke, spinal cord injuries, CVA, arthritic conditions or burns.

*Direction of strapping  
biomechanically controls the  
hand, wrist, and forearm.*

Left	Right	Size
NC25560	NC25561	X-Small
NC25562	NC25563	Small
NC25564	NC25565	Medium
NC25566	NC25567	Large

## One model, multiple applications

The versatile Ultra-Ortho™ eliminates the need of numerous other orthoses, including these Progress™ orthotics:



Progress™  
I.F. (Individual Finger)  
Ball Orthosis



Progress™  
Neutral Resting  
Orthosis



Progress™  
Neutral Hand Cone  
Orthosis



Progress™  
Dorsal Anti-Spasticity  
Orthosis



Progress™  
Palmar Bar Wrist  
Orthosis

To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call 800.821.9319




Your  
**COMPLETE SYSTEM**  
for **HAND** and  
Upper Extremity  
**THERAPY**

## TREAT **MULTIPLE** ISSUES.

Incorporating therapeutic  
movements of Flexion, Extension,  
Pronation and Supination,  
**FEPSim®** was developed to aid  
in the treatment of:

- Forearm, Wrist and Hand Fractures
- Neurological Conditions including Stroke
- Tendinopathies of the Upper Extremity
- Upper Extremity Weakness
- Osteoarthritis
- Complex Regional Pain Syndrome
- Pre and Post-Surgical Upper Extremity Conditions

## TWO **UNIQUE** MODELS. ONE **COMPLETE** SYSTEM.



To order, log in at **[www.ncmedical.com](http://www.ncmedical.com)** or call **800.821.9319**

# Big-Grip™ Telescoping Inspection Mirror

## Enhanced functionality for self-care.

- Four-way articulating mirror angles any direction for visibility.
- Telescoping stainless steel rod extends an additional 23".
- Built-up soft non-slip handle for easy grasping.
- Includes an adjustable safety lanyard.
- Mirror diameter: 2.75" (7cm).
- Product length ranges from 13" to 36" (33cm to 92cm).
- Requires (2) CR2032 batteries (included).
- Not made with natural rubber latex.

NC28743 Big-Grip™ Telescoping Inspection Mirror



*LED lights rotate 360° to illuminate dark areas*

To order, log in at **[www.ncmedical.com](http://www.ncmedical.com)** or call **800.821.9319**

Available in Virtual Showroom™

# Dressing Pal™

*Allows those with limited mobility to dress themselves more easily.*

- Combines a dressing stick and a shoehorn.
- Push-pull S-hook at one end, and a long-handled shoehorn at the opposite end.
- Hook end can be used to help put on clothing, or to reach clothes hangers on high closet rods.
- Other uses include retrieving clothing from floors, chairs or closets, and turning light switches on and off.
- Purple foam grip is soft and comfortable.



Dressing Pal™	Length
NC28581	24" (61 cm)
NC28582	30" (76 cm)



To order, log in at **[www.ncmedical.com](http://www.ncmedical.com)**  
or call **800.821.9319**

Available in Virtual Showroom™

# Shoehorns

Slide feet into shoes with comfort & ease.



## Big-Grip™ Shoehorns



- Built-up handle made of a latex-free\*, rubber-like material with flexible ribbing that adapts to any grip.
- Shaft is made of enameled steel.
- Textured surface prevents the foot from sticking to the shoehorn, helping it slide easily into the shoe.
- Handle diameter measures 1½" (3.8 cm).
- Available in three lengths.


\* Not made  
with natural  
rubber latex.



### Non-Slip Soft Grip

#### Big-Grip™ Shoehorns

NC65631	18" (46 cm)
NC65632	24" (61 cm)
NC65633	30" (76 cm)



## Easy-Slide Shoehorns

- Rubber-like handle for easy grasping.
- Enameled steel shaft.
- Available in three lengths.

#### Easy-Slide Shoehorns

NC28672	18" (46 cm)
NC28673	24" (61 cm)
NC28674	30" (76 cm)



## Foam Handle Shoehorns

- Built-up foam handle for gripping comfort.
- Enameled steel shaft.
- Available in three lengths.

#### Foam Handle Shoehorns

NC28608	18" (46 cm)
NC28609	24" (61 cm)
NC28610	30" (76 cm)

To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call **800.821.9319**

Available in Virtual Showroom™

**Norco®**  
Performance & Value

*Pill  
Pal™*

**SPLITTER**

***Safely split and  
easily store pills.***

- Split pills in half while keeping fingers safe.
- The top chamber includes a built-in safety shield, which hides the blade when the lid is opened.
- The bottom chamber can be used to store multiple pills and/or store unused portions of split pills.
- Pill splitter area measures 1.5" x 1.5" (3.8 x 3.8 cm).
- Device measures 3.5" x 2" (9 x 5 cm).
- Not made with natural rubber latex.

NC92101 Pill Pal™ Splitter



To order, log in at **[www.ncmedical.com](http://www.ncmedical.com)**  
or call **800.821.9319**  
Available in Virtual Showroom™

# Custom Hip Kits

## Custom Hip Kit #60

**Kit includes:**

- Norco® FeatherLite™ Reacher, 26" (66 cm)
- Norco® Molded Sock Aid with Two Handles
- Norco® Plastic Shoehorn with Hook
- Norco® Bath Sponge, Contoured

NC23000-60      Custom Hip Kit #60



## Custom Hip Kit #160

**Kit includes:**

- Easy-Slide Shoehorn, 18" (45 cm)
- Norco® FeatherLite™ II Reacher 26" (66 cm)
- Norco® Easy Pull Sock Aid
- After Total Hip Replacement Booklet

NC23000-160      Custom Hip Kit #160



## Custom Hip Kit #211

**Kit includes:**

- Norco® FeatherLite™ Reacher, 26" (66 cm)
- Norco® Molded Sock Aid with Two Handles
- Norco® Stainless Steel Shoehorn 24" (61 cm)
- Norco® Bath Sponge, Contoured
- Norco® Dressing Stick 27" (69 cm)
- Norco® Deluxe Elastic Shoe Laces, Black 27" (69 cm)

NC23000-211      Custom Hip Kit #211



*Visit our website for  
full hip kit product line.*

To order, log in at **www.ncmedical.com** or call **800.821.9319**

# Big-Grip™ Utensils

Cushioned grip keeps the utensil in the hand – even when wet.



**DISHWASHER  
SAFE  
STAINLESS  
STEEL**

**Non-Slip  
Soft Grip**



- Handle material is soft and latex-free\* with flexible ribbing that adapts to any grip.
- The special twist built into the metal shaft allows utensils to be bent to any angle.
- Weighted Utensil handles have 6 oz. (170 g) of added weight to help minimize tremors while eating.
- Coated teaspoon is helpful for those with biting reflexes or limited coordination.

## Big-Grip™ Utensils

NC65591	Fork
NC65595	Small Spoon
NC65593	Teaspoon, Uncoated
NC65597	Teaspoon, Coated
NC65594	Tablespoon
NC65590	Souper Spoon™
NC65592	Rocker Knife
NC65585	Serrated Rocker Knife
NC65599	Utensils Set (5)

## Big-Grip™ Weighted Utensils

NC65560	Fork
NC65561	Teaspoon
NC65562	Tablespoon
NC65564	Souper Spoon™
NC65565	Rocker Knife
NC65560-S	Utensils Set (5)



\*Not made with natural rubber latex.

To order, log in at **www.ncmedical.com**  
or call **800.821.9319**

Available in Virtual Showroom™

# Grippy™

NON-SLIP MATERIAL

## Get a Grip. Stop the Slip.



Trim to any  
size or shape.



*Round & Rectangle  
pre-cut mats  
COMING SOON!*



Not made  
with natural  
rubber latex.



- Flexible, non-slip material grips on both sides.
- Can be trimmed to any size or shape.
- Contains FDA approved antimicrobial additive.
- Reusable - easy to clean with soap and water.
- Easily removable from surfaces.
- Available in two widths: 8" (20cm) and 16" (41cm) — and two lengths: 2yd (1.8m) and 10yd (9.1m).

Grippy™	Roll Size
NC38800-1	8" x 2 yd (20cm x 1.8m)
NC38800-10	8" x 10 yd (20cm x 9.1m)
NC38801-1	16" x 2 yd (41cm x 1.8m)
NC38801-10	16" x 10 yd (41cm x 9.1m)

To order, log in at [www.ncmedical.com](http://www.ncmedical.com)  
or call **800.821.9319**

Available in Virtual Showroom™



**NEW COLORS**  
**GREAT QUALITY**  
**SAFETY-SEALED**  
**CONTAINERS**



**Rainbow™**

**EXERCISE PUTTY**

**Promote hand-muscle strength  
and increase endurance.**

- Exercising with Rainbow™ Exercise Putty helps promote increased range of motion, finger-flexion and tendon-gliding.
- Rainbow™ color-coded system associates resistance levels with color.
- The smaller portion cups are perfect for most small and medium sized hands.
- Composed of silicon rubber.



\* Not made with natural rubber latex.

**Latex Free\***

**Soft Orange**

NC98200-02	2 oz. (57 g)
NC98200-03	3 oz. (85 g)
NC98200-04	4 oz. (113 g)
NC98200-16	1 lb. (454 g)
NC98200-80	5 lbs. (2.3 kg)

**Medium-Soft Pink**

NC98210-02	2 oz. (57 g)
NC98210-03	3 oz. (85 g)
NC98210-04	4 oz. (113 g)
NC98210-16	1 lb. (454 g)
NC98210-80	5 lbs. (2.3 kg)

**Medium Green**

NC98220-02	2 oz. (57 g)
NC98220-03	3 oz. (85 g)
NC98220-04	4 oz. (113 g)
NC98220-16	1 lb. (454 g)
NC98220-80	5 lbs. (2.3 kg)

**Medium-Firm Blue**

NC98230-02	2 oz. (57 g)
NC98230-03	3 oz. (85 g)
NC98230-04	4 oz. (113 g)
NC98230-16	1 lb. (454 g)
NC98230-80	5 lbs. (2.3 kg)

**Firm Purple**

NC98240-02	2 oz. (57 g)
NC98240-03	3 oz. (85 g)
NC98240-04	4 oz. (113 g)
NC98240-16	1 lb. (454 g)
NC98240-80	5 lbs. (2.3 kg)

***Note:** The 3 oz. (90 cc of volume) is our most popular size Rainbow™ Exercise Putty and an ideal amount for large or small hands.*

## Putty Containers

- Additional empty containers for exercise putties.
- Cups can hold up to 5 oz. (141 g) of putty.
- Sold in a pack of 10 cups.



NC51025    **Putty Containers (Pack of 10)**

To order, log in at **www.ncmedical.com** or call **800.821.9319**

Available in Virtual Showroom™



# Air Putty®



Lightweight, velvety texture  
won't stick to skin or under fingernails.

- Perfect for gentle hand strengthening exercises.
- Unique composition of specialized silicon rubber is pleasant to feel.
- Maintains an even resistance while pinching and squeezing.
- Sold by volume, not by weight.
- Smaller portion volume is ideal for most small and medium sized hands.
- Not made with natural rubber latex.
- Free of the following contents: PABA, Fragrance, Animal protein, Gluten, Wheat, Shellfish, Tree nuts.

## Weight Conversion

60cc = 1.5oz  
90cc = 2.25oz  
120cc = 3oz  
400cc = 10oz (approx)  
1600cc = 2.5lbs (approx)

### Achieva™ Air Putty® Extra-Soft, Sky Blue

NC52405-60	60cc
NC52405-90	90cc
NC52405-120	120cc
NC52405-400	400cc
NC52405-800	1600cc

### Achieva™ Air Putty® Soft, Lavender

NC52410-60	60cc
NC52410-90	90cc
NC52410-120	120cc
NC52410-400	400cc
NC52410-800	1600cc

### Achieva™ Air Putty® Medium-Soft, Orchid

NC52415-60	60cc
NC52415-90	90cc
NC52415-120	120cc
NC52415-400	400cc
NC52415-800	1600cc

### Achieva™ Air Putty® Medium, Pink

NC52420-60	60cc
NC52420-90	90cc
NC52420-120	120cc
NC52420-400	400cc
NC52420-800	1600cc

To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call **800.821.9319**

Available in Virtual Showroom™

# Rainbow™

LATEX FREE  
EXERCISE BAND

A premium  
**Latex Free\* Exercise Band**  
without a premium price!

LEVEL 1

LEVEL 2

LEVEL 3

LEVEL 4

LEVEL 5



Color-Coded Resistance Levels with Functional Colors™

- No-Slip Grip • No Powder • No Odor



**Rainbow™ Latex-Free\* Exercise Band**

**Professional**

**50 yd Box**

One 50 yd (45 m) roll

NC95401-50	Level 1 - Light	Orange
NC95402-50	Level 2 - Medium	Pink
NC95403-50	Level 3 - Heavy	Green
NC95404-50	Level 4 - Extra Heavy	Blue
NC95405-50	Level 5 - Super Heavy	Purple

**Retail**

**Multipacks - Three 5 ft (1.5 m) lengths of band per pack**

NC95410	Light Resistance Pack	Orange, Pink, Green	Levels 1, 2 and 3
NC95411	Heavy Resistance Pack	Green, Blue, Purple	Levels 3, 4 and 5

**Complete Multipack - Five 5 ft (1.5 m) lengths of band per pack**

NC95400	Orange, Pink, Green, Blue, Purple	Levels 1, 2, 3, 4 and 5
---------	-----------------------------------	-------------------------

**Latex Free\***

\*Not made with natural rubber latex.

To order, log in at  
**www.ncmedical.com** or call **800.821.9319**

Available in Virtual Showroom™

# Outriggers

## Phoenix<sup>®</sup> & Base 2<sup>™</sup>

**Back in Stock!**

- For fabrication of dynamic and static progressive splints.
- Available in convenient kits or separate components to create customized contoured outriggers.
- Lightweight aluminum components can be bent by hand or with pliers.
- Wing nuts, thumb screws and line guides made from nylon and stainless steel.
- Finger slings are made of ultra-suede and are available pre-tied or with grommets.



### Phoenix<sup>®</sup> and Base 2<sup>™</sup> Outrigger Kits

		PCS
NC12820	Base 2 Single Finger Extension Kit	7
NC12821	Base 2 Kit #2 Multi Digit Extension Kit	14
NC12840	Base 2 Outrigger Starter Kit	63
NC12760	Phoenix Outrigger Kit	15
NC12763	Phoenix Extended Outrigger Kit	12

### Base 2<sup>™</sup> Outrigger Aluminum Components

		Length	QTY
NC12830-3	Straight Aluminum	2" (5 cm)	2
NC12830-1	Straight Aluminum	5" (13 cm)	3
NC12830-0	Straight Aluminum	7" (18 cm)	3
NC12830-2	Curved Multi-digit	4" (10 cm)	2
NC12830-10	Curved Outrigger	5" (13 cm)	3
NC12830-6	Curved Outrigger	7" (18 cm)	3
NC12830-8	Curved Outrigger w/foot	5" (13 cm)	3
NC12830-7	Curved Outrigger w/foot	7" (18 cm)	3

### Outrigger Accessoriess

		Length	QTY
NC12338	Finger Slings w/ Grommets	2" (5 cm)	10
NC12397	Finger Slings w/ Grommets	5" (8 cm)	10
NC12410	Finger Slings w/ Grommets	3.5" (9 cm)	10
NC12398	Finger Slings w/ Grommets	4" (10 cm)	10
NC12399	Ultra-Suede Pre-tied Sling	3" (8 cm)	10
NC12401	Ultra-Suede Pre-tied Sling	4" (10 cm)	10
NC12830-4	Base 2 Accessories Kit		13
NC12762	Rubber Band Posts (Thumb screws)		10
NC12830-9	Wing Nuts		5
NC12830-11	All Plastic Line Guides		10
NC12494	Outrigger Line Connectors		100
NC12760-1	Phoenix Slotted Pulleys w/hex wrench		5
NC12760-2	Radial Outrigger Bar		10



To order, log in at [www.ncmedical.com](http://www.ncmedical.com) or call 800.821.9319



# Virtual Showroom

E-COMMERCE AFFILIATE PROGRAM

Join the Affiliate Program. Get an instant **online store**.



**Virtual Showroom™** is an easy to install, fully-functional e-commerce platform designed specifically for rehab healthcare providers. Virtual Showroom™ integrates seamlessly into your existing website and gives your clients direct access to the products their therapists have recommended. It is simple to create your own Virtual Showroom™ and it requires no programming skills.

## *Highlights:*

- Clients stay on your website when shopping for recommended products.
- Customizable: Select Featured Products, Categories, and color-match to your website.
  - Adapts to mobile devices and tablets.

## *Benefits:*

- No annual, monthly, or transaction fees.
  - No inventory required.
- Up to 25% commissions on all product sales.
  - Competitive prices for all products.
- North Coast handles all customer service inquiries.
- Full product warranties and our 100% Satisfaction & Return Policy for all sales.

Sign-up at **[www.mynpcmstore.com/virtual-showroom](http://www.mynpcmstore.com/virtual-showroom)**  
or call **800.821.9319**

# Rx Cart<sup>TM</sup> 2.0 *As easy as writing a prescription.*

RxCart<sup>TM</sup> 2.0 allows the therapist to specify the right products, with the proper size, weight, resistance level, and brand. Your clients no longer need to wander aimlessly through a confusing website hoping to get the right products with proper fitment. This puts the expert decisions in the hands of you, the expert.


Harness the power of the therapist recommendation.

Simple to use on a tablet or smartphone while consulting with patients.

Makes it easy for your patient to get the specific therapy products you recommend.

## **Features include:**

- Analytics to track sales to the location and the therapist.
- Save the product bundles you recommend most.
- View status, send reminders for ignored RxCarts.
- Organizations can proactively track and influence therapist's behavior, and reward those who are using the tools best.
- Up to 25% commissions on all product sales.



**You specify the products. Your patient clicks a link in their email.  
In just a few easy steps, patients purchase your list of recommended products.**

**Bring the power of Virtual Showroom<sup>TM</sup> into your daily practice.**

Sign-up at [www.myncmstore.com/rx-cart](http://www.myncmstore.com/rx-cart)  
or call **800.821.9319**



## Our Company

For 50 years, North Coast Medical has maintained an outstanding reputation in the global rehabilitation market. We are more than just a wholesale distributor and manufacturer of products, equipment and supplies. North Coast is a strategic supply chain partner, providing innovative opportunities for our customers and business partners.

Now a second-generation family business, we continue to invest in our future while steadfastly adhering to strong ethics and integrity. We invite you to browse our website and discover what is *more than expected* at North Coast Medical.

Thousands of medical and rehabilitative products are available at [www.ncmedical.com](http://www.ncmedical.com) – log in to view your professional pricing. Orders can also be placed by phone at **800.821.9319**, or by fax at **877.213.9300**.

## Proprietary Brands

Throughout the years, we've developed and/or acquired medical devices with trademarks and brand names including:



# North Coast

MEDICAL & REHABILITATION PRODUCTS

*More Than Expected™*