



BLUE DIAMOND™
DRY NEEDLES



YOUR **ONE** BRAND FOR
DRY NEEDLING!



YOUR **ONE** BRAND FOR
DRY NEEDLING!

BLUE DIAMOND™
DRY NEEDLES



8 Needle Lengths / 8 Needle Diameters / 20 Needle Combinations

- Features an oversized, conductive aluminum etched handle, ensuring superior control when utilizing enhanced target-specific techniques.
- Japanese high-grade surgical steel needle body. Designed to allow maximum flexion and deeper penetration.
- Medical-grade silicone lubricated needle body engineered to glide effortlessly and reduce friction – remarkably gentle for improved patient comfort.
- Each box contains 20 sleeves of 5 needles (100 pieces per box).
- FDA cleared.


*Easy to peel
medical dialysis
paper backing*



To order, log in at **www.ncmedical.com** or call **800.821.9319**

*Oversized, conductive, etched
aluminum handle*

*Double-clamped tube for easy
one-handed release*



BLUE DIAMOND™
DRY NEEDLES

NEEDLE DIAMETER COLOR CHART



0.16mm
Red



0.18mm
Coral



0.20mm
Blue



0.22mm
Pink



0.25mm
Purple



0.30mm
Brown



0.40mm
Tan



0.45mm
Orange

0.16 to 0.22mm Diameter

	Diameter	Length
NCBD-1613	0.16mm	13mm
NCBD-1813	0.18mm	13mm
NCBD-2013	0.20mm	13mm
NCBD-2213	0.22mm	13mm

0.25mm Diameter

	Diameter	Length
NCBD-2513	0.25mm	13mm
NCBD-2530	0.25mm	30mm
NCBD-2540	0.25mm	40mm
NCBD-2550	0.25mm	50mm
NCBD-2560	0.25mm	60mm
NCBD-2575	0.25mm	75mm
NCBD-25100	0.25mm	100mm
NCBD-25125	0.25mm	125mm

0.30mm Diameter

	Diameter	Length
NCBD-3030	0.30mm	30mm
NCBD-3040	0.30mm	40mm
NCBD-3050	0.30mm	50mm
NCBD-3060	0.30mm	60mm
NCBD-3075	0.30mm	75mm
NCBD-30100	0.30mm	100mm

0.40 to 0.45mm Diameter

	Diameter	Length
NCBD-40100	0.40mm	100mm
NCBD-45125	0.45mm	125mm

*Reserved for licensed medical
professionals and approved
resellers only.*



Rainbow™

LATEX FREE
EXERCISE TUBING

**Consistent Resistance
Poundage Force**
while stretching!



Norco®

Rainbow™ Latex-Free* Exercise Tubing

Professional

100 ft Box

One 100 ft (30m) roll

NC95421-100	Level 1 - Light	Orange
NC95422-100	Level 2 - Medium	Pink
NC95423-100	Level 3 - Heavy	Green
NC95424-100	Level 4 - Extra Heavy	Blue
NC95425-100	Level 5 - Super Heavy	Purple

Retail

Multipacks

Three 5 ft (1.5m) lengths of tubing per pack

NC95430	Light Resistance Pack	Orange, Pink, Green	Levels 1, 2 and 3
NC95431	Heavy Resistance Pack	Green, Blue, Purple	Levels 3, 4 and 5

Color-Coded Resistance Levels with Functional Colors™

- No Residue
- No Powder
- No Odor

Latex Free*

*Not made with natural rubber latex.

To order, log in at
www.ncmedical.com or call **800.821.9319**

Available in Virtual Showroom™

Rainbow™

LATEX FREE
EXERCISE BAND

A premium
Latex Free* Exercise Band
without a premium price!

LEVEL 1

LEVEL 2

LEVEL 3

LEVEL 4

LEVEL 5



Color-Coded Resistance Levels with Functional Colors™

- No-Slip Grip • No Powder • No Odor



Rainbow™ Latex-Free* Exercise Band

Professional

50 yd Box

One 50 yd (45 m) roll

NC95401-50	Level 1 - Light	Orange
NC95402-50	Level 2 - Medium	Pink
NC95403-50	Level 3 - Heavy	Green
NC95404-50	Level 4 - Extra Heavy	Blue
NC95405-50	Level 5 - Super Heavy	Purple

Retail

Multipacks - Three 5 ft (1.5 m) lengths of band per pack

NC95410	Light Resistance Pack	Orange, Pink, Green	Levels 1, 2 and 3
NC95411	Heavy Resistance Pack	Green, Blue, Purple	Levels 3, 4 and 5

Complete Multipack - Five 5 ft (1.5 m) lengths of band per pack

NC95400	Orange, Pink, Green, Blue, Purple	Levels 1, 2, 3, 4 and 5
---------	-----------------------------------	-------------------------

Latex Free*

*Not made with natural rubber latex.



To order, log in at
www.ncmedical.com or call **800.821.9319**

Available in Virtual Showroom™



Rainbow™

Exercise Band Mini Loops

Short, continuous band loops are ideal for upper and lower body exercises.


- Rainbow™ Exercise Band Mini Loops are seam-free and flat.
- Mini Loops measure 2" (5 cm) x 11" (28 cm).
- Available in five color-coded resistances.
- Contains latex.
- Sold individually.



Rainbow™ Exercise Band Mini Loops

	Resistance	Color
--	------------	-------

NC91241-01	Level 1	Orange
NC91242-01	Level 2	Pink
NC91243-01	Level 3	Green
NC91244-01	Level 4	Blue
NC91245-01	Level 5	Purple



To order, log in at www.ncmedical.com
or call **800.821.9319**

Available in Virtual Showroom™

Exercise Balls

Economical burst-resistant balls provide versatile exercise workouts.

- Tested to support up to 550 lbs. (250 kg).
- Helps to improve core stabilization by training the spinal column muscles as the body works to stay in balance.
- Balls are available individually in retail box with handle.
- Not intended for heavy lifting.
- Not made with natural rubber latex.

Ideal for clinic or home use.

Exercise Balls

	Color	Diameter
NC50100	Blue	45 cm (17¾")
NC50101	Red	55 cm (21½")
NC50102	Teal	65 cm (25½")
NC50103	Dark Blue	75 cm (29½")
NC50104	Silver	85 cm (33½")




Mini Ball

Strengthen key muscle groups, core health, endurance, posture and balance.

- High quality, 9" (23 cm) diameter therapy ball intended for exercising and strengthening the abdomen, back, hips, knees, shoulders, and neck.
- Anti-burst vinyl material deflates slowly if punctured.
- Designed to remain elastic and squishy when fully inflated.
- Supports up to 330 lbs. (150 kg).
- Not made with natural rubber latex.

Includes exercise poster.



NC50109 Mini Ball

To order, log in at www.ncmedical.com or call 800.821.9319

Available in Virtual Showroom™


Economy Shoulder Pulley

Highest quality in the market!

- Increase shoulder range of motion in the clinic or at home.
- Large, easy-grip handles.
- Non-slip web door strap.
- Simple to adjust cord length.
- Smooth action pulley mechanism.



Economy
Shoulder Pulley
NC52080



Universal Door Bracket

Easy on, easy off
over-the-door bracket.

- Bracket can remain in place on the top of the door, ready for the next exercise.
- Made of smooth, coated metal to help prevent door scratches.
- Bracket width is 1³/₄" (4 cm), which fits over commercial and most standard doors.
- Includes two EVA foam pads which help secure the bracket around doors with narrow widths.



Use self-adhesive
foam pads for
narrower
door widths.



1-3/4"
(4.4 cm)



NC52081 Universal Door Bracket

To order, log in at www.ncmedical.com or call 800.821.9319

Available in Virtual Showroom™


Shoulder Pulley

Everything needed to perform
range of motion exercises at home.

- Ideal for home use, to independently improve range of shoulder movement and coordination following shoulder injury or surgery.
- Can be used for exercises while seated, standing, or lying down.
- Sets up in minutes by positioning the door strap over the top along the hinged side of any standard door.
- Pre-assembled and ready to use.
- The door strap fits snugly between the door and frame and will not scratch surfaces.
- For optimal attachment to door, use the Norco® Universal Door Bracket (sold separately).
- Not made with natural rubber latex.



*Door strap fits snug
between door
and frame.*



NC52069 Shoulder Pulley



Retail ready box includes:

- Exercise booklet with step-by-step instructions and photos
- Door strap
- 6½ ft. (2.4 m) cord
- Pulley mechanism and comfortable, ergonomically designed handles

To order, log in at **www.ncmedical.com** or call **800.821.9319**

Available in Virtual Showroom™

Stretch Strap

Safe and effective deep,
gradual stretching.



ISOTONICS
ISOMETRICS
PROLONGED STRETCH
INCREASED R.O.M.
FLEXIBILITY

- Designed with ten multi-positioned grip loops for use in deep, gradual stretching of the major muscle groups.
- Ideal for individuals seeking unassisted or assisted stretching exercises to relieve sore or achy muscles.
- Stretch more safely and effectively than is possible unaided.
- Grip loops can be used to set stretch improvement goals.
- Made of non-stretch webbing.
- Machine washable.
- Not made with natural rubber latex.

*Retail ready
packaging includes
an instruction booklet
for quick reference.*



Stretch Strap

	Model	Color	Length
NC84527	Standard	Purple	6' 4" (1.9 m) for people under 6' (1.8 m)
NC84528	Long	Blue	7' 5" (2.3 m) for people over 6' (1.8 m)

To order, log in at **www.ncmedical.com** or call **800.821.9319**

Available in Virtual Showroom™

Common Trigger Points



BODYtool™

TRIGGER POINT SELF MASSAGER

Repair Yourself

- Designed to target trigger points and muscle knots in the body to assist with temporary release of muscular tension.
- **Textured grip handles** allow better control for effectively targeting trigger points.
- Features **seven large knobs**, strategically placed to control pressure to the desired areas for trigger point release and breaking up of myofascial tissue.

NC69700 Body Tool™

**Body Tool™
Retail Stand
NC69700-RS**



Retail display rack holds 12 Body Tool™ massagers.



Easy-to-follow instruction guide and poster included!

To order, log in at **www.ncmedical.com** or call **800.821.9319**

Available in Virtual Showroom™




EVA Foam Roller

Professional quality foam roller designed for heavy and repeated use in multi-user environments.

- Great for proprioceptive balance and stability exercise or to roll out trigger points, knots and areas of increased muscle density.
- Use for strengthening, flexibility and stretching exercises.
- Provides comfortable support — gentle on tight muscles while lying supine, yet functional enough for light massage rolling.
- Closed-cell, textured, high-density EVA foam technology resists breaking down and allows for easy cleaning between clients.
- Built with high-grade materials for long-lasting performance.
- Roll measures 36" (91 cm) long with a 6" (15 cm) diameter.

EVA Foam Roller
NC84800-01



Not made with
natural rubber latex.

To order, log in at www.ncmedical.com or call 800.821.9319

Available in Virtual Showroom™

Lumbar Roll


Provides comfortable support to help improve posture and relieve lower back pain.

- Adjustable, non-elastic strap allows users to secure the Lumbar Roll to standard chair backs and car seats at the desired height.
- Innovative strapping design allows the Lumbar Roll to maintain position in the hollow of the back, or be repositioned for maximum comfort.
- Made of polyurethane foam with a removable polyester cover that is easy to clean.
- The removable cover includes a zipper guard for easy opening and closing of the zipper.
- Available in **Standard** and **Firm** densities.
- Rolls measure 11" (28 cm) long with a 5" (13 cm) diameter.
- Not made with natural rubber latex.



Lumbar Roll

	Foam Density
NC89095	Standard
NC89096	Firm



Strap "tunnel" allows the strap to be tightened firmly to maintain correct positioning.



Zipper guard for snag-free zipper use.



SCAN
to learn more!

To order, log in at **www.ncmedical.com** or call **800.821.9319**

Available in Virtual Showroom™



MUSCLE AIDTAPE™

ARM YOUR MUSCLE

BEAT THE PAIN

- Relieve muscle pain, soreness and cramping.
- Provide support for joints and muscles.
- Reduce swelling, inflammation and bruising.
- Increase blood circulation/flow.
- Water-resistant and breathable.
- Not made with natural rubber latex.



Available in:

Single rolls: 2" x 16' (5 cm x 5 m)

Bulk rolls: 2" x 105' (5 cm x 32 m)

Display Cartons	Color	Quantity
MT1642-MVBG	Beige	6 Rolls
MT1642-JDBK	Black	6 Rolls
MT1642-RPOR	Orange	6 Rolls
MT1642-MVPK	Pink	6 Rolls
MT1642-JDNB	Navy Blue	6 Rolls
MT1642-JDBL	Blue	6 Rolls
MT1642-RPUS	USA	6 Rolls

Bulk Rolls	Color
MT1052-BG-BULK	Beige
MT1052-BK-BULK	Black
MT1052-BL-BULK	Blue



Also Available in
10" Long Precut Strips

Strengthen Performance.
Relieve Muscle Pain.
Facilitate Rapid Recovery.



Precut strips: 2" x 10" (5 cm x 25.5 cm)

Precuts (20)	Color	Bulk Precuts (128)	Color
MTPCBG-20	Beige	MTPCBG-128	Beige
MTPCBK-20	Black	MTPCBK-128	Black
MTPCBL-20	Blue	MTPCUS-128	USA
MTPCUS-20	USA		

**MUSCLE
AIDTAPE™**

Longer lasting, skin-friendly support.



Reid Priddy
Olympic Gold Medalist




To order, log in at **www.ncmedical.com** or call **800.821.9319**

Available in Virtual Showroom™

RELIEVE PAIN NATURALLY

Premium Hemp Products

- Helps relieve pain safely and quickly, reducing inflammation that may be caused by injury and overuse of muscles/joints.
- Pain relief formulated by Ayurveda PH.D. doctors and hemp experts.
- Organic, USA-grown hemp, all natural plant-based formula contains Ayurveda herbs and medical grade essential oils.
- Produced in an FDA, USDA Organic, ISO, and WHO-GMP certified facility to ensure the highest quality and standards.



Full-Spectrum
Hemp Extract
Pain Lotion
1,000 mg
NC70724
50.3 mL / 1.7 oz
Pump




Full-Spectrum
Hemp Extract
Pain Oil
500 mg
NC70722
10 mL / 0.3 fl oz
Roll-on



Broad-Spectrum
(THC-Free)
Hemp Extract Pain Oil
1,000 mg
NC70723
10 mL / 0.3 fl oz
Roll-on

To order, log in at **www.ncmedical.com**
or call **800.821.9319**
Available in Virtual Showroom™



CLINICAL

MAXIMUM PERFORMANCE PAIN RELIEF.

- Provides fast, long-lasting relief from sprains, strains, bruises, arthritis, muscle aches, joint and back pain.
- Specially formulated penetration enhancers to provide fast pain relief.
- Contains 10% menthol, dimethyl sulfone (MSM) and glucosamine.
- Sweat-resistant.
- Paraben-free, no dyes or fragrances.
- Non-greasy and non-staining.



Stopain® Clinical - Singles Qty

NC80091-02	Roll-On	3 fl. oz. (88 ml)	1
NC80091-03	Spray	4 fl. oz. (118 ml)	1
NC80091-01	Gel Tube	4 fl. oz. (118 ml)	1

Stopain® Clinical - Larger Size Gel Pump Qty

NC80091-05	32 fl. oz. (946 ml)	1
NC80091-06	1 gallon (3.78 liters)	1
NC80091-04	16 fl. oz. (473 ml)	1

Stopain® Clinical - Starter Pack Qty

NC80091-07	Trial Packets w/ Display Box	(100) 5 gram packets
------------	------------------------------	----------------------



To order, log in at www.ncmedical.com or call 800.821.9319

Available in Virtual Showroom™



Eco-Stim™ TENS + EMS

Combo unit provides safe, drug-free pain relief.

- TENS, EMS and Massage wave forms.
- 22 treatment modes. 6 modes allow Time, Frequency and Pulse Width customization.
- Large, easy-to-read backlit display.
- 40 levels of output intensity.
- Safety feature locks intensity level.
- Dual channel 120mA (at 500 Ohm load).
- Output frequency:
2 - 100 Hz
Pulse width:
100- 300 μ s per pulse
- Includes:
(4) 2"x 2" square reusable electrodes,
(2) lead wires,
(3) AAA batteries,
(1) belt clip, and
(1) instruction manual.

NC89481
Eco-Stim™ TENS + EMS

*Use belt clip for
on-the-go pain relief.*



Not made with
natural rubber latex.

*Includes 48 page
instruction manual.*



To order, log in at **www.ncmedical.com** or call **800.821.9319**

Available in Virtual Showroom™



Latex Free*

Eco-Stim™

Neurostimulation Electrodes



- Premium quality Carbon/Hydrogel for maximum comfort.
- Reusable up to 20 times when stored as recommended.
- Compatible with all NMS and TENS devices.
- Personalized packaging for patient data usage.
- FDA cleared for OTC purchase.
- 2mm DIN pigtail connector.
- Available in a variety of shapes and sizes in both foam or cloth
- Not made with natural rubber latex.

Cloth



Foam



*Personalized packaging
for patient data usage.*

*4 electrodes
per pack*




**2" x 2"
Square**

Foam: NC89292-1
Cloth: NC89293-1




**2" x 4"
Rectangle**

Foam: NC89292-3
Cloth: NC89293-3



**2" diameter
Round**

Foam: NC89292-2
Cloth: NC89293-2



**2 3/4" diameter
Round**

Cloth: NC89293-4

To order, log in at
www.ncmedical.com
or call **800.821.9319**
Available in Virtual Showroom™



Polyvinyl Cold Packs

**Professional cold packs remain flexible
and comfortable down to -4°F (-20°C).**



NC70595
Oversize • 22 x 13 in. (56 x 33 cm)



NC70596
Neck Contour • 23 x 8 in. (58 x 20 cm)



NC70594
Standard • 15 x 12 in. (38 x 30 cm)



NC70593
Half Size
9½ x 5½ in. (24 x 14 cm)

- Designed for professional or at-home use.
- Effective relief of pain, discomfort or swelling after injury for up to 30 minutes.
- Filled with a premium, non-toxic gel that remains pliable and withstands freezing solid.
- Guaranteed not to split, crack or leak.
- One year manufacturer warranty.
- Not made with natural rubber latex.



To order, log in at
www.ncmedical.com
or call **800.821.9319**
Available in Virtual Showroom™



Polyurethane Cold Packs

Premium cold packs conform to body contours at freezing temperatures.

- Designed to withstand heavy clinic use.
- Effective relief of pain, discomfort or swelling after injury for up to 30 minutes.
- Filled with a premium, eco-friendly gel that is pliable at freezing temperatures.
- Guaranteed not to split, crack or leak.
- One year manufacturer warranty.
- Not made with natural rubber latex.



NC70598

Oversize • 22 x 13 in. (56 x 33 cm)



NC70597

Standard • 15 x 12 in. (38 x 30 cm)



NC70599

Neck Contour • 23 x 8 in. (58 x 20 cm)



To order, log in at www.ncmedical.com or call **800.821.9319**

Available in Virtual Showroom™

Wireless "On-the-Go" Technology

Intelligent iontophoresis.



Technologically advanced,
self-contained iontophoresis
drug delivery system.

- Fully self-contained.
- Intelligent biofeedback microprocessor.
- Hydrogel dispersive technology—**no saline required.**
- Patented "shut-off" technology.
- SmartPower™ LED.
- Detects and measures skin impedance to optimize drug delivery.

Active area: 11.4 cm²
Dosage: 40 or 80 mA*min
Skin Adhesive: Hypoallergenic
Maximum Buffer Capacity: Ag/AgCl Stabilizer
Fill Volume: 2.0cc
Smart Power: 6V with intelligent microprocessor
Packaging: 6 ActivaPatch® Kits, Alcohol Wipes,
Patient Information Cards
Treatment time: 2.5 Hours

ActivaPatch® IntelliDose™ 2.5
NC89188



**80 mA*min dose
in only 2.5 hours.**

* Study data on file

no saline needed.

To order, log in at www.ncmedical.com or call 800.821.9319

ActivaPatch® IontoGo™ technology.

No external power source to connect, no leadwires,
no dose controller needed, power shut down feature and
no saline required.

ACTIVApatch®
IontoGo™ 4.0

Extended wear patch with an average drug delivery dose of 80 mA*min in approximately 4 hours.



Active area: 11.4 cm²
Dosage: 80 mA*min
Skin Adhesive: Hypoallergenic
Fill Volume: 2.0cc
Packaging: 6 IontoGo™ 4.0,
6 Alcohol Wipes
Treatment Time: 4 Hours

**NO
SALINE
REQUIRED**

**ActivaPatch®
IontoGo™ 4.0
NC89189**

**FULLY
SELF-
CONTAINED**



ACTIVApatch®
IontoGo™ 12.0

Extended wear patch with an average drug delivery dose of 80 mA*min in approximately 12 hours.



Active area: 11.4 cm²
Dosage: 80 mA*min
Skin Adhesive: Hypoallergenic
Fill Volume: 2.0cc
Packaging: 6 IontoGo™ 12.0,
6 Alcohol Wipes
Treatment Time: 12 Hours

**NO
SALINE
REQUIRED**

**ActivaPatch®
IontoGo™ 12.0
NC89191**

**FULLY
SELF-
CONTAINED**



**HOW SMART
IS YOUR
WIRELESS IONTO?**

Smarter Iontophoresis™

* Study data on file

To order, log in at **www.ncmedical.com** or call **800.821.9319**

Wired Technology




ADVANTAGES

- Conforms extremely well to uneven treatment sites.
- Maximum strength, precise buffering with new Carbon-Ag/AgCl technology.
- Large ground electrode surface area lowers sensation and increases comfort.
- Highly absorbant drug matrix.
- Large surface adhesion area.



TRIVARION® IS AVAILABLE IN THE FOLLOWING SIZES:



ADVANTAGES

- Two controls operate dose and current, with optional push-knob or turn capabilities.
- High-contrast display indicates dose, time and current.
- More power than dual-channel devices for better penetration; up to 80 volts.
- Automatic current ramp-up and ramp-down for maximum comfort. Dose range from 0 to 80 mA-min with 0.1 mA-min increments.
- Time calculations are performed automatically and recalculated if dose and current settings are changed, or with pause/restart.



ActivaDose® II Controller
NC89186

To order, log in at www.ncmedical.com or call 800.821.9319



Buffered Iontophoretic Delivery Kit



Each electrode has a large surface adhesion area.

Variety of sizes and shapes offer great flexibility for iontophoresis treatments.

- Made of non-sensitizing, hypoallergenic, medical grade materials.
- Delivers ions with a pH balance up to 80 mA-min.
- Maximum strength, precise buffering with unique carbon- Ag/AgCl technology.
- The highly absorbent drug matrix has no messy hydrogel.
- Large ground electrode surface area lowers sensation and increases comfort.
- Each kit contains 12 individual treatment sessions that include: one activity electrode; one return electrode; one alcohol prep pad.

Flexible. Thin. Contouring. Comfortable.



Competitor's product is not as thin or flexible.



Advanced, technology-based manufacturing for a thinner, more pliable electrode.



Easy-to-open packaging.

Not made with natural rubber latex.



Small
Square 1.5 cc
NC89250B



Medium
Square 2.5 cc
NC89251B



Large
Square 4.0 cc
NC89252B



Butterfly 2.0 cc
NC89255B

To order, log in at www.ncmedical.com or call **800.821.9319**

Goniometers

**Accurate, economical
and easy-to-read.**

International Goniometer 360°

- Linear measurements up to 7" and up to 18.5 cm.
- 3 rings for measuring ROM:
Outer ring 0-90° in 10° increments bilaterally
Middle ring 0-180° in 1° increments bilaterally
Inner ring 0-360° in 20° increments.

NC70102

International Goniometer 180°



- Linear measurements up to 12.5" and up to 32 cm.
- Outer ring 0-180° in 1° increments bilaterally.

NC70119

180° Plastic Goniometer

- Linear measurements up to 4" and up to 11 cm.
- 180° scale is marked in 5° increments.

NC70100



Devore® Pocket Goniometer

- Measures either direction from 0-180° in easy-to-read 1° increments.

NC70115

To order, log in at
www.ncmedical.com
or call **800.821.9319**

Body Weight Scale

Digital scale accurately measures body weight up to 440 lbs.

- Easily switch between pounds (lbs) or kilograms (kg).
- Large digital numbers are easy to read.
- Provides an audible beep when weight is confirmed.
- Data displays for 10 seconds after stepping off scale.


Large, easy-to-read digital display.

- Low battery/ overload indicator.
- Tempered safety glass platform.
- Dimensions: 13.5" x 12.25" x 1" (34.3 x 31 x 2.5 cm).
- 3 AAA batteries are included.

Step on scale to turn on. Scale shuts off automatically.



Sleek, low-profile design adds an elegant look to clinic or home.



Retail-ready packaging.

NC66860 Body Weight Scale

To order, log in at www.ncmedical.com or call 800.821.9319

Available in Virtual Showroom™

Gait Belts

A safe, easy way to transfer clients.



Rainbow
Pastel



Yellow

- Gait Belts can be used as walking and transferring belts.
- Available in four colors:
White, Blue, Rainbow Pastel, or Yellow.
- **Yellow** gait belt indicates high-risk of falls due to decreased safety awareness.
- Made of soft, heavy-duty webbing that is easy to clean.
- Belts measure 2" (5.1 cm) wide.
- Machine washable.
- Measure waist to determine size.



Blue

*Teeth lock
for a secure hold.*

Gait Belts

	Color	Length	Fits Waists
NC91370	Yellow	54" (137 cm)	15" to 51" (38 to 130 cm)
NC91380	White	54" (137 cm)	15" to 51" (38 to 130 cm)
NC91393	Blue	54" (137 cm)	15" to 51" (38 to 130 cm)
NC91391	Rainbow Pastel	54" (137 cm)	15" to 51" (38 to 130 cm)
NC91371	Yellow	60" (152 cm)	15" to 54" (38 to 137 cm)
NC91386	White	60" (152 cm)	15" to 54" (38 to 137 cm)
NC91394	Blue	60" (152 cm)	15" to 54" (38 to 137 cm)
NC91381	Rainbow Pastel	60" (152 cm)	15" to 54" (38 to 137 cm)
NC91373	Yellow	72" (183 cm)	15" to 68" (38 to 173 cm)
NC91387	White	72" (183 cm)	15" to 68" (38 to 173 cm)
NC91383	Blue	72" (183 cm)	15" to 68" (38 to 173 cm)
NC91392	Rainbow Pastel	72" (183 cm)	15" to 68" (38 to 173 cm)



White

To order, log in at www.ncmedical.com or call **800.821.9319**

Polyvinyl Gait Belt

**Easy-to-clean belt helps
prevent cross-contamination.**



- Made of vinyl-coated webbing that can be easily cleaned.
- Features a metal buckle that provides a firm grip to prevent slippage.
- Polyvinyl material is sonic-welded and sealed, eliminating crevices for germs to accumulate.
- Belts measure 2" (5.1 cm) wide.
- Available in three lengths.

Polyvinyl Gait Belts
Length

NC91355	54" (137cm)
NC91356	60" (152cm)
NC91357	72" (183cm)



Real-world motion capture *without the cost.*

Capture. Analyze. Document. Implement.

PROPRIETARY TECHNOLOGY

12 years in development, 5 engineering teams, 60 engineers, and over \$10 million invested. The result is the unique patented technology of Digitsole®. The Digitsole® **Mov-Scan™ pods** efficiently capture gait and lower extremity biomechanical movement with remarkable accuracy, validity and reproducibility.



Digitsole® **Mov-Scan™ pods capture over 100,000 points-of-data in under 1 minute** and offer the opportunity to quickly assess and evaluate both quantitative and qualitative movement. Digitsole® allows for real-time, real-world practical solutions for your clients.

HARDWARE

Each Digitsole® Pro Assessment Kit is portable and includes a full set of Digitsole® Pro insoles and a pair of Mov-Scan™ pods for capturing and recording steps, strides and foot orientation while walking, running or jumping.

Digitsole's proprietary algorithms calculate lower extremity, spatiotemporal kinematic and advanced parameters with remarkable precision and accuracy. Mov-Scan™ pods utilize the world's smallest combination of a 3-D Accelerometer, Gyroscope, and Inertial Measurement Unit (IMU).




SOFTWARE

Digitsole® software is intuitive and extremely easy to use. Data is organized, easily accessible and presented to empower the healthcare practitioner in order to quickly analyze, document, and implement therapeutic solutions.

Included in every Digitsole® system is the ability to determine which Stride AI™ Foot Orthotic, model, type and postings are most appropriate for your client.





*Analyze and follow
your patients' activity in
real-world environments.*



CAPTURE DATA FOR ASSESSMENT IN 1-3 MINUTES

The algorithms capture over 32 parameters in real-world conditions. This solution provides an objective interpretation to improve the clinical assessment of patients.



In-depth walking analysis 1 minute



In-depth running analysis 2 minutes



Performance Rehab Jump Module 3 minutes

Digitsole Pro's **Mov-Scan™** monitors sports injury rehabilitation, improves athletic performance, and helps prevent injury by assessing risks and tracking results.



Video Module

Simultaneously capture video and movement data.

Digitsole® Pro Analytical System

Walk & Run Modules, Hardware Kit, Site License,
2 Users, insole & posting recommendations

NC51551

Digitsole® Basic Analytical System

Walk & Run insole recommendation, Hardware
Kit, Site License, 2 Users, insole & posting
recommendations

NC51551-01

Digitsole® Analytical Modules:

Digitsole® Walk Module

NC51552-04

Digitsole® Walk & Run Module

NC51552-05

Performance Rehab Jump Module

NC51552-02

Video Module

NC51552-03

Digitsole® Add-ons:

Additional Hardware Kit

NC51550

2 Additional Users Pack

NC51552-06

5 Additional Users Pack

NC51552-01



SCAN
to learn more!

To order, log in at **www.ncmedical.com** or call **800.821.9319**

Orthotics that support *the way you move.*

Recommended by Digitsole® Pro advanced analytics.



HEEL STRIKE

RUN

Stride AI[™] Run, Heel Strike

Deep heel cup with maximum cushion and control.

Run Heel Strike Part #	Size	QTY	Men's Shoe Size	Women's Shoe Size
NC51561-XS	XS	1 pair	3½ to 5	4½ to 6
NC51561-S	S	1 pair	5½ to 7	6½ to 8
NC51561-M	M	1 pair	7½ to 9	8½ to 10
NC51561-L	L	1 pair	9½ to 11	10½ to 12
NC51561-XL	XL	1 pair	11½ to 13	12½ to 14
NC51561-XXL	XXL	1 pair	13½ to 15	14½ to 15

ALL STRIDE AI[™] RUN MODELS FEATURE:

- ✓ HEAT-MOLDABLE ECO FOAM BASE
- ✓ CONTOURED BASE TO ENHANCE YOUR PERFORMANCE
- ✓ CUSTOMIZABLE POSTING ZONES

Stride AI[™] Run, Midfoot

Unique mid foot support with cushion.

Run Midfoot Part #	Size	QTY	Men's Shoe Size	Women's Shoe Size
NC51562-XS	XS	1 pair	3½ to 5	4½ to 6
NC51562-S	S	1 pair	5½ to 7	6½ to 8
NC51562-M	M	1 pair	7½ to 9	8½ to 10
NC51562-L	L	1 pair	9½ to 11	10½ to 12
NC51562-XL	XL	1 pair	11½ to 13	12½ to 14
NC51562-XXL	XXL	1 pair	13½ to 15	14½ to 15



MIDFOOT

FOREFOOT

Stride AI[™] Run, Forefoot

Unique forefoot shape to maximize forefoot propulsion.

Run Forefoot Part #	Size	QTY	Men's Shoe Size	Women's Shoe Size
NC51560-XS	XS	1 pair	3½ to 5	4½ to 6
NC51560-S	S	1 pair	5½ to 7	6½ to 8
NC51560-M	M	1 pair	7½ to 9	8½ to 10
NC51560-L	L	1 pair	9½ to 11	10½ to 12
NC51560-XL	XL	1 pair	11½ to 13	12½ to 14
NC51560-XXL	XXL	1 pair	13½ to 15	14½ to 15



POSTINGS/WEDGES

See website for details.

To order, log in at www.ncmedical.com

or call **800.821.9319**

Available in Virtual Showroom[™]

Orthotics that support *the way you move.*

Recommended by Digitsole® Pro advanced analytics.



COMFORT

ALL STRIDE AI™ WALK MODELS FEATURE:

- ✓ HEAT-MOLDABLE ECO FOAM BASE
- ✓ HIGH PERFORMANCE TOP COVER
- ✓ CUSTOMIZABLE POSTING ZONES



SUPPORT

PERFORMANCE



POSTINGS/WEDGES

See website for details.

WALK

Stride AI™ Walk, Comfort

Low firmness for kids, elderly, lightweight patients, or those who need minimal biomechanical support.

Walk Comfort Part #	Size	QTY	Men's Shoe Size	Women's Shoe Size
NC51563-XS	XS	1 pair	3½ to 5	4½ to 6
NC51563-S	S	1 pair	5½ to 7	6½ to 8
NC51563-M	M	1 pair	7½ to 9	8½ to 10
NC51563-L	L	1 pair	9½ to 11	10½ to 12
NC51563-XL	XL	1 pair	11½ to 13	12½ to 14

Stride AI™ Walk, Support

Medium firmness for moderate biomechanical support.

Walk Support Part #	Size	QTY	Men's Shoe Size	Women's Shoe Size
NC51565-XS	XS	1 pair	3½ to 5	4½ to 6
NC51565-S	S	1 pair	5½ to 7	6½ to 8
NC51565-M	M	1 pair	7½ to 9	8½ to 10
NC51565-L	L	1 pair	9½ to 11	10½ to 12
NC51565-XL	XL	1 pair	11½ to 13	12½ to 14

Stride AI™ Walk, Performance

Reinforced shell for maximum support during activity.

Walk Performance Part #	Size	QTY	Men's Shoe Size	Women's Shoe Size
NC51564-XS	XS	1 pair	3½ to 5	4½ to 6
NC51564-S	S	1 pair	5½ to 7	6½ to 8
NC51564-M	M	1 pair	7½ to 9	8½ to 10
NC51564-L	L	1 pair	9½ to 11	10½ to 12
NC51564-XL	XL	1 pair	11½ to 13	12½ to 14

To order, log in at **www.ncmedical.com**
or call **800.821.9319**

Available in Virtual Showroom™



Virtual Showroom

E-COMMERCE AFFILIATE PROGRAM

Join the Affiliate Program. Get an instant **online store**.



Virtual Showroom™ is an easy to install, fully-functional e-commerce platform designed specifically for rehab healthcare providers. Virtual Showroom™ integrates seamlessly into your existing website and gives your clients direct access to the products their therapists have recommended. It is simple to create your own Virtual Showroom™ and it requires no programming skills.

Highlights:

- Clients stay on your website when shopping for recommended products.
- Customizable: Select Featured Products, Categories, and color-match to your website.
 - Adapts to mobile devices and tablets.

Benefits:

- No annual, monthly, or transaction fees.
 - No inventory required.
- Up to 25% commissions on all product sales.
 - Competitive prices for all products.
- North Coast handles all customer service inquiries.
- Full product warranties and our 100% Satisfaction & Return Policy for all sales.

Sign-up at **www.mynpcmstore.com/virtual-showroom**
or call **800.821.9319**

Rx CartTM 2.0 *As easy as writing a prescription.*

RxCartTM 2.0 allows the therapist to specify the right products, with the proper size, weight, resistance level, and brand. Your clients no longer need to wander aimlessly through a confusing website hoping to get the right products with proper fitment. This puts the expert decisions in the hands of you, the expert.


Harness the power of the therapist recommendation.

Simple to use on a tablet or smartphone while consulting with patients.

Makes it easy for your patient to get the specific therapy products you recommend.

Features include:

- Analytics to track sales to the location and the therapist.
- Save the product bundles you recommend most.
- View status, send reminders for ignored RxCarts.
- Organizations can proactively track and influence therapist's behavior, and reward those who are using the tools best.
- Up to 25% commissions on all product sales.



You specify the products. Your patient clicks a link in their email.
In just a few easy steps, patients purchase your list of recommended products.

Bring the power of Virtual ShowroomTM into your daily practice.

Sign-up at www.myncmstore.com/rx-cart
or call 800.821.9319



Our Company

For 50 years, North Coast Medical has maintained an outstanding reputation in the global rehabilitation market. We are more than just a wholesale distributor and manufacturer of products, equipment and supplies. North Coast is a strategic supply chain partner, providing innovative opportunities for our customers and business partners.

Now a second-generation family business, we continue to invest in our future while steadfastly adhering to strong ethics and integrity. We invite you to browse our website and discover what is *more than expected* at North Coast Medical.

Thousands of medical and rehabilitative products are available at www.ncmedical.com – log in to view your professional pricing. Orders can also be placed by phone at **800.821.9319**, or by fax at **877.213.9300**.

Proprietary Brands

Throughout the years, we've developed and/or acquired medical devices with trademarks and brand names including:



North Coast

MEDICAL & REHABILITATION PRODUCTS

More Than Expected™