

Eco-Stim™ TENS + EMS

22
Treatment
Modes

40
Levels of
Output

Combo unit provides safe, drug-free pain relief.


- TENS, EMS and Massage wave forms.
- 22 treatment modes. 6 modes allow Time, Frequency and Pulse Width customization.
- Large, easy-to-read backlit display.
- 40 levels of output intensity with safety lock feature.
- Dual channel 120mA (at 500 Ohm load).
- **Output frequency:** 2 - 100 Hz
- **Pulse width:** 100- 300 μ s per pulse
- Not made with natural rubber latex.
- Includes:
 - (4) 2" x 2" square reusable electrodes,
 - (2) lead wires,
 - (3) AAA batteries,
 - (1) belt clip, and
 - (1) instruction manual.



NC89481
Eco-Stim™ TENS + EMS



Use belt clip for on-the-go pain relief.



Includes 48 page instruction manual.

For more information, contact your North Coast Medical representative or visit www.ncmedical.com