

BAILEY

Model 809 Training Stairs Our Smallest Stairs - 30" Wide - Pass Thru

Staircase Features:

- 30" wide access, with 30" wide x 24" deep Platform at 24"H.
- 500 LB capacity, with solid hardwood and engineered hardwood construction.
- Single handrail design at 28" height.
- Step height 6" (2 each), and 4" (3 each), tread depth 10", with anti-slip tread and platform surfaces
- Includes removable Bus Step for public transportation simulation
- Overall Footprint is 55"L x 33"W
- Max height at upper platform handrail is 40".
- Minimum required ceiling height 89" to 105".
- Ships "knocked down", (197 LBS) assembly required

