

ActivaPatch® Performance Features

Wireless “On-the-Go” Technology

No external power source to connect, no lead wires, no dose controller needed.



Features	ActivaPatch® IntelliDose™ 2.5	ActivaPatch® IontoGo™ 4.0	ActivaPatch® IontoGo™ 12.0
Treatment Time	2.5 hours	4 hours	12 hours
Dosage	40 or 80 mA*min	80 mA*min	80 mA*min
Drug Reservoir Fill Volume	2.0cc	2.0cc	2.0cc
Hydrogel Dispersive Technology No Saline Required	✓	✓	✓
Pull-Tab Activation	✓	✓	✓
Patented “shut-off” technology	✓	✓	✓
Intelligent biofeedback microprocessor	✓		
SmartPower™ LED	✓		