

# BAILEY

## Model 3100 Balance Platform

### ***Balance Platform Features:***

- 36" x 36" Carpeted Platform.
- Platform height 2" including carpet.
- 500 LB capacity, with solid wood platform, and steel posts and handrails.
- Dual handrails adjust from 26" to 44" high.
- 1 1/2" Diameter Handrails, with Plunger pull and release locking height adjustment
- 30" between uprights allows manual and some power wheelchair access.
- Ships "knocked down", assembly required
- Ships via UPS, weight 89 lbs.

**Model 3100**  
*36" x 36" Platform*  
*Covered with HD Low Pile Carpet*



**Model 3100 Shown With Model 725**  
*Multiple Configurations Possible*  
*Use as mini Training Steps*

