

# YOUR ONE BRAND FOR DRY NEEDLING!

- Designed with an oversized, etched aluminum handle for superior control.
- Japanese high-grade surgical steel needle body.
- Medical-grade silicone lubricated needle for improved patient comfort.
- Each box contains 20 sleeves of 5 needles (100 pieces per box).



**8 Needle Lengths**  
**8 Needle Diameters**  
**20 Needle Combinations**



## BLUE DIAMOND™ DRY NEEDLES

### 0.16 to 0.22mm Diameter

|           | Diameter | Length |
|-----------|----------|--------|
| NCBD-1613 | 0.16mm   | 13mm   |
| NCBD-1813 | 0.18mm   | 13mm   |
| NCBD-2013 | 0.20mm   | 13mm   |
| NCBD-2213 | 0.22mm   | 13mm   |

### 0.25mm Diameter

|            | Diameter | Length |
|------------|----------|--------|
| NCBD-2513  | 0.25mm   | 13mm   |
| NCBD-2530  | 0.25mm   | 30mm   |
| NCBD-2540  | 0.25mm   | 40mm   |
| NCBD-2550  | 0.25mm   | 50mm   |
| NCBD-2560  | 0.25mm   | 60mm   |
| NCBD-2575  | 0.25mm   | 75mm   |
| NCBD-25100 | 0.25mm   | 100mm  |
| NCBD-25125 | 0.25mm   | 125mm  |


### 0.30mm Diameter

|            | Diameter | Length |
|------------|----------|--------|
| NCBD-3030  | 0.30mm   | 30mm   |
| NCBD-3040  | 0.30mm   | 40mm   |
| NCBD-3050  | 0.30mm   | 50mm   |
| NCBD-3060  | 0.30mm   | 60mm   |
| NCBD-3075  | 0.30mm   | 75mm   |
| NCBD-30100 | 0.30mm   | 100mm  |

### 0.40 to 0.45mm Diameter

|            | Diameter | Length |
|------------|----------|--------|
| NCBD-40100 | 0.40mm   | 100mm  |
| NCBD-45125 | 0.45mm   | 125mm  |

### NEEDLE DIAMETER COLOR CHART



*Oversized, conductive, etched aluminum handle*



*Double-clamped tube for easy one-handed release*



Distributed by: