

Bariatric Treatment Table



MODEL 6190

- Plywood top and solid hardwood legs
- 6-leg design with triple bolted corner legs
- Built-in middle leg levelers
- Steel reinforced frame
- Laminate full length shelf with H-braces
- Pillow included
- 2" firm foam padding
- 600 lbs. capacity under normal use

SPECIFICATIONS

6190

78" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Options



077- Positioning Strap

034- Paper Dispenser & Holder Combo

040- Paper Cutter

030- Paper Dispenser

Finishes



Natural



Dark Cherry

Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black