



# MULTI-PURPOSE TABLES

Our extremely versatile units can be used for a variety of functions:



**MODEL 6293**  
W/ DRAWER AND CASTERS

- **Hand Therapy**  
Can be used by a therapist and patient for hand therapy.
- **Mobile Therapy Table**  
13" height range allows the patient to use while sitting, standing, ambulating, or in a wheelchair.
- **Mobile Computer Table**  
Lightweight, easily maneuverable for use anywhere.

## MODELS:

- 6292** Mobile Table
- 6293** Mobile Table w/Drawer
- 6294** Mobile Table w/Cut Out

MODELS	LENGTH	WIDTH	HEIGHT	w/CASTERS	DRAWER	FRONT CUT OUT
6292	37"	19"	29"- 42"	30"- 43"	N/A	N/A
6293	37"	19"	29"- 42"	30"- 43"	16"W x 17"D	N/A
6294	37"	24"	29"- 42"	30"- 43"	N/A	19"W x 6"D

### Models 6292, 6293, 6294 Multi-Purpose Tables



- Dual gas spring combined with an unweighted table top allows the table top to be effortlessly lifted to required height.
- Dual hand control levers on each side releases and lowers the height to desired level.
- Brushed aluminum dual telescoping legs with black powder coated U-shaped base.
- Thick, High Pressure Laminate tops (1"), finished in Folkstone Gray with black T-molding edge banding and rounded corners.
- 2" diameter twin-wheel black nylon locking casters and 5/8"-1" leveler glides are included with each unit for mobile or stationery configurations.
- 250 lbs static weight capacity.



**MODEL 6294**  
W/ FRONT CUT OUT AND CASTERS