

HOOK SPLINT.®

A sleek, low profile Relative Motion Orthosis.
Relative motion made extremely simple.



Secure a custom fit with adjustable locking tabs.



Relative motion is a technique which places one or more fingers in a relatively extended or flexed position compared to the others.

Allows the entire hand to stay functional while achieving treatment objectives.

Designed as a ready-to-use alternative to thermoplastics.

Formed by a metallic band with a smooth silicone finish, and a locking tab on each end.

Waterproof, UV coating protects the outer surface and ensures long-term patient comfort.

Indications include:

Extensor tendon/sagittal band repair,
PIP extension or flexion contracture,
boutonniere deformity, trigger finger,
Dupuytren's disease, tendon adhesions,
digit pain/stiffness.

	Size	Length	Width
Singles			
NC58774-1	SM/MED	11" (28 cm)	.39" (1 cm)
NC58775-1	LG	11.8" (30 cm)	.47" (1.2 cm)
Case of 10			
NC58774-10	SM/MED	11" (28 cm)	.39" (1 cm)
NC58775-10	LG	11.8" (30 cm)	.47" (1.2 cm)

SCAN for VIDEO



Distributed by:

