



## Treatment Tables Specifications

### Table Item #30841

Table Description: 2-Section  
Overall Dimension: 27 in. W, 81 in. L  
Cushion Dimension: 27 in. W, 77 in. L  
Head Section: 27" x 27", raises 90 degrees or lowers 45 Degrees  
Foot Section: 27" x 48"  
Lower Position From Ground: 19"  
Upper Position From Ground: 36"  
Total Weight: 221.34 lbs



### Table Item #30805

Table Description: 7-Section  
Overall Dimension: 27 in. W, 85 in. L  
Cushion Dimension: 27 in. W, 75 in. L  
Head Section: 15" x 13", raises 30 degrees or lowers 30 degrees  
Armrests: Range of 6 in. vertical adjustment  
Center Section: 27" x 26", will raise 15 degrees  
Foot Section: 27" x 33", raises 90 degrees  
Lower Position From Ground: 21"  
Upper Position From Ground: 38"  
Total Weight: 285.50 lbs



### Table Item #30364

Table Description: Traction Table  
Overall Dimension: 27 in. W, 82 in. L  
Cushion Dimension: 27 in. W, 71 in. L  
Head Section: 27" x 15", raises 90 degrees or lowers 45 degrees  
2 Center Section: 27" x 12" (two parts) – Center section slide 6 ½" open  
Foot Section: 27" x 28", raises 70 degrees  
Lower Position From Ground: 19"  
Upper Position From Ground: 36"  
Total Weight: 288.14 lbs



### **Table Item #32060**

Table Description: Bo Bath  
Overall Dimension: 40 in. W, 83 in. L  
Cushion Dimension: 40 in. W, 83 in. L  
Head Section: 40" x 26", raises 90 degrees or lowers 25 degrees  
Foot Section: 40" x 56"  
Lower Position From Ground: 18"  
Upper Position From Ground: 36"  
Total Weight: 321.43 lbs



### **Table Item #32089**

Table Description: 90 Degree Table  
Overall Dimension: 27 in. W, 84 in. L  
Cushion Dimension: 27 in. W, 78 in. L  
Head Section: 27" x 14", raises 75 degrees or lowers 90 degrees  
Center Section: 27" x 24"  
Foot Section: 27" x 36", raises 90 degrees  
Lower Position From Ground: 18"  
Upper Position From Ground: 36"  
Total Weight: 221.34 lbs



The Overall Dimension is the maximum space the table will occupy in any configuration (up, down, or in-between). This can be from foot pedal-to-foot pedal or cushion-to-cushion, or anywhere in between.