


Big-Grip™ Telescoping Inspection Mirror

Enhanced functionality for self-care.

- Four-way articulating mirror angles any direction for visibility.
- Telescoping stainless steel rod extends an additional 23".
- Built-up soft non-slip handle for easy grasping.
- Includes an adjustable safety lanyard.
- Mirror diameter: 2.75" (7cm).
- Product length ranges from 13" to 36" (33 cm to 92 cm).
- Requires (2) CR2032 batteries (included).
- Not made with natural rubber latex.

NC28743 Big-Grip™ Telescoping Inspection Mirror



LED lights rotate 360° to illuminate dark areas

Distributed by: