

# Norco® Finger Extension Tube

## Reduce flexion contractures.

- The shape and design of the Norco® Finger Extension Tube provides a low stress, high duration capability, ideal for flexion contractures.
- Circumferential pressure helps restrict fluid to aid in edema reduction.
- Tip of the finger should be visible. If needed, trim either end.
- Orthotic should fit comfortably, a firm fit is not required for effective results.
- Not recommended for flexion contractures greater than 50°.
- Made of dense neoprene, polyester, and poly/cotton tape.
- Not made with natural rubber latex.



	Size	PIP Joint Circumference
NC79068-1	X-Small	1-3/4" to 2" (4.5 to 5.1cm)
NC79068-2	Small	2" to 2-1/2" (5.1 to 6.4cm)
NC79068-3	Medium	2-1/2" to 2-3/4" (6.4 to 7.0cm)
NC79068-4	Large	2-3/4" to 3" (7.0 to 7.6cm)
NC79068-5	X-Large	3" to 3-1/2" (7.6 to 8.9cm)



*Poly/cotton tape provides comfort and durability.*

