




GLADIATOR™

Designed in conjunction with Hand, OT and PT professionals to produce a **firm** orthosis that provides **active range of motion**. The **formfitting** orthosis can be used when treating wrist sprains, TFCC injuries, and postsurgical distal radius fractures.

Short Wrist Orthosis

Breathable fabric



3 removable metal stays

Biomechanical Advantage

- Counterforce straps
- Formfitting
- 1/2 forearm length
- Removable stays
- Breathable
- Lightweight

Comfort Cool® Gladiator™ Short Wrist Orthosis

Part #	Size	Hand	Wrist Circumference
NC86640	Small	Left	6" to 7" (15 to 18 cm)
NC86641	Small	Right	6" to 7" (15 to 18 cm)
NC86642	Medium	Left	7" to 8" (18 to 20 cm)
NC86643	Medium	Right	7" to 8" (18 to 20 cm)
NC86644	Large	Left	8" to 9" (20 to 23 cm)
NC86645	Large	Right	8" to 9" (20 to 23 cm)
NC86646	X-Large	Left	9" to 10" (23 to 25 cm)
NC86647	X-Large	Right	9" to 10" (23 to 25 cm)

Counterforce straps



"The Gladiator Short Orthosis is a game changer in the off-the-shelf market, offering versatile protection and support for various wrist conditions, including post-surgical and chronic arthritis, as well as cumulative trauma disorders like carpal tunnel syndrome, while being easily adjustable to meet individual patient needs and maximize occupational performance."

– Greg Pitts, OTD, OTR/L, CHT



Visit www.ncmedical.com or contact your North Coast Medical sales representative for details. Toll Free: **800.821.9319**