



INTELECT® LEGEND 2 COMBO POWERED BY VMS

Electrotherapy Simplified

The Intellect® Legend 2 is a two-channel or four-channel combo electrotherapy system used with or without an optional Cart. Intellect Legend 2 features Chattanooga's proprietary VMS waveform which enables stronger, more efficient muscle contractions¹ while the suggested protocols reduce treatment setup time to help improve clinic workflow.

KEY FEATURES

- VMS waveform achieves a stronger muscle contraction while increasing patient comfort and decreasing fatigue.^{2,3,4}
- Suggested protocol settings (SPS). Select parameters by selecting indications instead of waveforms.
- Save up to 12 protocol shortcuts on the home screen custom protocols. Begin treatments within 2 button presses.
- Hi resolution 7" touch screen
- On board memory and software upgrade feature
- 12 waveforms for electrical stimulation: VMS, Interferential, Premodulation, Asymmetric Biphasic, Symmetric Biphasic, High Volt, Microcurrent, Russian, VMS Burst, HAN, Direct Current, VMS FR
- 1 MHz and 3 MHz ultrasound




Variable Muscle Stimulation (VMS)

INTELECT® LEGEND 2 COMBO STANDARD ACCESSORIES:

PART NUMBER	DESCRIPTION
81-5010	Intellect Legend 2 Two Channel Combo Set
81-5011	Intellect Legend 2 Four Channel Combo Set
79967	Carbon Electrodes
42198	Dura-Stick+ 5cm Square Electrodes QTY 4
70010	Stim Lead Wires CH 1&2
70011	Stim Lead Wires CH 3&4 (4CH combo only)
12-10648	Nylatex Strap
15-1062	5cm ² Ultrasound Applicator
14679	Power Cord

INTELECT® LEGEND 2 OPTIONAL ACCESSORIES:

15-1136	Intellect Legend 2 Cart
14-1086	Intellect Legend 2 Battery
15-0160	1cm ² applicator
15-0161	2cm ² applicator
15-0163	10cm ² applicator
4248-12	Ultrasound Gel (12 Bottles)
79967	6x8 cm Carbon Electrodes (QTY 4)
15-1140	USB drive



[chattanoogarehab.com/us/clinical-electrotherapy](http://chattanooga rehab.com/us/clinical-electrotherapy)

chattanooga®

CLINICAL ELECTROTHERAPY DEVICE COMPARISON



	Vectra Neo	Intellect Legend 2	Intellect Transport 2
Channels	2 or 4	2 or 4	2
Warranty	2 Years	2 Years	2 Years
Save Custom Protocols	200	25	15
Modular Design	4 Modules		
Remote Control	✓		
Low Level Laser	✓		
sEMG	✓		
Color Screen	✓	✓	✓
Includes VMS	✓	✓	✓
Touch Screen	✓	✓	
Clinical Protocols	✓	✓	
Magnetic Cart Attachment		✓	✓
Waveforms	12	12	5

enovis™

T 800.321.9549 F 800.936.6569

DJO, LLC
5919 Sea Otter Place, Suite 200 | Carlsbad, CA 92010 | U.S.A.
enovis.com/chattanooga

Copyright © 2023 DJO, LLC
MKT00-13033 Rev B

1. Bellew, Allen, Biefnes, Grantham, Miglin, and Swartzell, 2018

2. Adams, Harris, Woodard, and Dudley, 1993; Binder-Macleod, Halden, and Jungles, 1995; Gorgey and Dudley, 2008

3. Bowman and Baker, 1985

4. Cayot, and Bellew. J Sports Rehabil. 2022

Individual results may vary. Neither DJO, LLC nor any of the Enovis companies dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.