

# Dema Grip® Gentle

## Tubular support with light compression

- Designed to provide approximately half the amount of pressure compared to regular Dema Grip® Compression Stockinette.
- Use single layered, or doubled over on itself with a 4/5" to 1 1/5" (2 to 3cm) exterior layer edge overlap.
- Can be used over wounds to secure bandages, or under orthoses for added protection, compression or padding.
- Made of blended cotton and spandex.
- Size based on flat bandage width.
- Rolls measure 33' (10m) long.
- Not made with natural rubber latex.



**NOW AVAILABLE IN  
BLACK!**



*Provides a lighter  
amount of compression  
for sensitive skin types.*

### Dema Grip® Gentle Compression Stockinette

Natural	Black	Size	Width	Length
NC42700	NC42800	A	2" (5.1 cm)	33' (10 m)
NC42701	NC42801	B	2 1/2" (6.4 cm)	33' (10 m)
NC42702	NC42802	C	2 3/4" (7 cm)	33' (10 m)
NC42703	NC42803	D	3" (7.6 cm)	33' (10 m)
NC42704	NC42804	E	3 1/2" (8.9 cm)	33' (10 m)
NC42705	NC42805	F	4" (10 cm)	33' (10 m)
NC42707	NC42807	G	4 3/4" (12 cm)	33' (10 m)
NC42708		J	6 3/4" (17 cm)	33' (10 m)
NC42709		K	8" (20 cm)	33' (10 m)
NC42710		L	10" (25 cm)	33' (10 m)
NC42711		M	12 3/4" (32 cm)	33' (10 m)