



# Pill Pal™

## SPLITTER

*Safely split & easily store pills.*

- Split pills in half while keeping fingers safe.
- The top chamber includes a built-in safety shield, which hides the blade when the lid is opened.
- The bottom chamber can be used to store multiple pills and/or store unused portions of split pills.
- Pill splitter area measures 1.5" x 1.5" (3.8 x 3.8 cm).
- Device measures 3.5" x 2" (9 x 5 cm).
- Not made with natural rubber latex.

Blister  
Card  
retail  
package



NC92101 Pill Pal™ Splitter



Visit [www.ncmedical.com](http://www.ncmedical.com) or contact your North Coast Medical sales representative for details.  
Toll-free: 800.821.9319