

# Outriggers

## Phoenix<sup>®</sup> & Base 2<sup>™</sup>

### Back in Stock!

- For fabrication of dynamic and static progressive splints.
- Available in convenient kits or separate components to create customized contoured outriggers.
- Lightweight aluminum components can be bent by hand or with pliers.
- Wing nuts, thumb screws and line guides made from nylon and stainless steel.
- Finger slings are made of ultra-suede and are available pre-tied or with grommets.



Phoenix <sup>®</sup> and Base 2 <sup>™</sup> Outrigger Kits			PCS
NC12820	Base 2 Single Finger Extension Kit		7
NC12821	Base 2 Multi-Digit MP Extension Kit		14
NC12840	Base 2 Outrigger Starter Kit		63
NC12760	Phoenix Outrigger Kit		15
NC12763	Phoenix Extended Outrigger Kit		12

Base 2 <sup>™</sup> Outrigger Aluminum Components		Length	QTY
NC12830-3	Straight Aluminum	2" (5 cm)	2
NC12830-1	Straight Aluminum	5" (13 cm)	3
NC12830-0	Straight Aluminum	7" (18 cm)	3
NC12830-2	Curved Multi-digit	4" (10 cm)	2
NC12830-10	Curved Outrigger	5" (13 cm)	3
NC12830-6	Curved Outrigger	7" (18 cm)	3
NC12830-8	Curved Outrigger w/foot	5" (13 cm)	3
NC12830-7	Curved Outrigger w/foot	7" (18 cm)	3

Outrigger Accessories		Length	QTY
NC12338	Finger Slings w/ Grommets	2" (5 cm)	10
NC12397	Finger Slings w/ Grommets	5" (8 cm)	10
NC12410	Finger Slings w/ Grommets	3.5" (9 cm)	10
NC12398	Finger Slings w/ Grommets	4" (10 cm)	10
NC12399	Ultra-Suede Pre-tied Sling	3" (8 cm)	10
NC12401	Ultra-Suede Pre-tied Sling	4" (10 cm)	10
NC12830-4	Base 2 Accessory Kit		13
NC12762	Rubber Band Posts (Thumb screws)		10
NC12830-9	Wing Nuts		5
NC12830-11	Line Guides		10
NC12494	Outrigger Line Connectors		100
NC12760-1	Phoenix Slotted Pulleys w/hex wrench		5
NC12760-2	Phoenix Radial Outrigger Bar		1



Distributed by: