

The Ideal Needle for Dry Needling Just Got Better

As a Dry Needling practitioner, you know it requires deep tissue work to stimulate myofascial trigger points to ease pain, cure injuries and speed-up recovery. It requires a lot of dexterity and deeper needle insertions. For years practitioners improvised with needles that were not made for Dry Needling. Red Coral worked closely with educators and practitioners to create the first purpose-built dry needle. We continue to work with educators and practitioners to innovate and evolve our needles. Myotech 2.0 Elite Dry Needles[®] are the latest evolution that features two more dry needle innovations.

Thicker, conductive handle with etched grip


- + Better dexterity when handling, particularly on deeper insertions
- + Quicker to e-stim to confirm placement and accelerate healing

Stronger German surgical steel that stays sharp

- + Won't bend or break during deeper insertions
- + Allows for multiple insertions with each needle

Full needle parylene coating for green chemistry lubrication

- + Easier to penetrate patient skin and tissue for tougher trigger points and deeper insertions
- + A near painless patient experience



Longer Needle Lengths

- + Allows for deeper trigger point penetrations

Holding-tab and glue free guide tube

- + One-handed needle insertions
- + Reduced waste and time

Thicker Needle Diameters

- + Easier to handle
- + Won't break on deeper insertions
- + Allows for longer needle lengths and deeper trigger point penetrations

Micro-channeling

- + Easier to penetrate patient skin and tissue for tougher trigger points and deeper insertions
- + A near painless patient experience

MYOTECH[®] DRY NEEDLES



“Having practiced and taught Myofascial Dry Needling since 1986, I have used a great variety of needles, without doubt Myotech is the best I have used.”

Andrew Gallagher,
Consultant Physiotherapist
(Pain Management), Senior
Lecturer at The Myofascial
Pain Study Centre,
Melbourne, Australia

Better needles. Better needling treatments.
Better patient and practitioner experiences.[®]



Even Better Needles for Better Patient and Practitioner Experiences

The first generation Myotech needle featured several innovations that made a marked difference in patient comfort and practitioner ease of use during dry needling treatments.

- > German surgical steel that does not break or bend during deeper insertions, and stays sharp through multiple insertions
- > Thicker etched handle and bigger needle sizes that make the needles easier to handle and go deeper
- > Full needle parylene coating that allows for near painless insertions
- > Microscopic grooves along the needle shaft that provide micro-channels to reduce surface area. Significantly reducing the force needed for insertions
- > Holding tab free, glued guide tube for easy one-handed aim, snap, and tap insertions

The 2nd generation Myotech needle adds two new, patent pending technologies and features that further limit patient discomfort and enhance practitioner efficiency and effectiveness.

- > A double clamped, 30% thicker handle that improves dexterity and strength
- > Glue-free guide that uses a dimple inside the tube for more efficient aim and tap insertions

MYOTECH[®]
DRY NEEDLES

Myotech 2.0 Elite Dry Needles[®]

The first purpose-built dry needle that was designed and made for the strength and dexterity needed for deeper insertions while limiting patient discomfort and improving practitioner efficiency. The Myotech 2.0 now has an even thicker handle and features a glue-free guide tube.

Features

- + Stronger German surgical steel
- + Needle surface micro-channeling
- + Full needle parylene coating
- + Thick conductive handle with etched grip
- + High-visibility, holding tab and glue free guide tube

Available Needle Gauges and Lengths

GAUGE	LENGTHS										
	13mm	15mm	30mm	40mm	50mm	60mm	70mm	75mm	100mm	135mm	
.16											
.18			✓	✓							
.20		✓	✓	✓							
.22		✓									
.25			✓	✓	✓						
.30			✓	✓	✓	✓		✓	✓		
.35					✓	✓		✓			
.40									✓		
.45								✓		✓*	

* Available with tab-free, glued clear guide tube only



Contact your representative today:
(800) 821-9319 | www.ncmedical.com

