



Prescribed for Progress™



ISO7000R Recumbent Bike Product Specifications

Order Code: ISO7011R-INT

FEATURES

Bi-directional movement and resistance	Pedaling in forward and reverse provides versatility, helps prolong exercise and promotes reciprocal muscle balance
Low starting resistance	6 watts
Resistance range	191 levels of resistance - 20 levels adjustable in .1 increments
Resistance system	3 phase combination generator & eddy current brake - largest watt range in the industry
Large, self-righting pedals	Strap adjusts easily to fit the user
Handlebars	Handlebars on seat and console provide a variety of hand grip options
Easy to transport	Integrated transport wheels
Heart rate monitoring	Contact and telemetric (chest strap required)

ACCESSIBILITY & ADJUSTABILITY



Step-through seating	Adjustable step-through allows easy entry and safe and easy access for all users
Direct wheelchair access	Seat slides off for direct wheelchair access - ADA compliant
Premium seat	Seat adjusts up/down, forward/back, swivels and reclines
360° Therapist Assist Pedals	Adjust the seat forward and back from any point around the bike
Seat Recline	Recline to any point within a 23 degree range
Seat Swivel	360 degree swivel that locks at each 90 degrees for safety
Seat height adjustment	Select seat height anywhere from 18.5"-23" (46cm-58cm)
Seat removal handle	Integrated handle on seat back makes it easy to remove the seat
Accessible adjustment levers	Reach all adjustment levers from both sides of the seat
Zero-sloping forward/back seat adjustment	Adjust forward and back on level monorail for optimum control and safe positioning

INTELLI-FIT™ CONSOLE

7" LCD full color touch screen display	Touch screen with clear, simple screens and large, readable font.
Audible beep	Yes
Tactile buttons	Tactile buttons in addition to touch screen
Resistance levels	191 levels of resistance - 20 levels adjustable in .1 increments
Display feedback	Time, level, distance, METs, heart rate, calories, RPM, watts
One touch quick start	Yes
Control METs - Constant Work Program	Set constant METs level and the workload remains constant throughout the program.
Iso-Strength Program	Isokinetic program for safe, accommodating resistance
Power Burst	8-second burst of high intensity training
Cool down	Automatic cool down begins after each exercise is completed
Administrative mode	User defined settings, metric or US units of measure and other custom settings
User-defined console languages	9 languages - English, Spanish, French, German, Russian, Japanese, Chinese, Italian, Dutch
Download results to USB flash drive	Plug USB flash drive into the Fit-Key® slot to save workout results in a .txt file

INTELLI-FIT™ CONSOLE (cont.)

CSAFE port and 8-volt DC power port	Yes			
USB charger	For phones and other electronic devices			
Water bottle holder, accessory tray, reading rack	Integrated			
Personal Cooling Fan	Integrated dual fans with three speeds			
Programs	Quick Start	Constant Work - METs	Hills	Heart Fit Test & Training
	Manual	Constant Work - Watts	Fit-Quik®	Power Fit Test & Training
	Heart Rate	Iso-Strength	Random	Stress Test

TECHNICAL SPECIFICATIONS

Power Requirements	Self-generating; low-voltage AC adapter, included
Maximum user weight capacity	450 lbs. (204kg)
Unit weight	243 lbs. (106kg)
Unit size	58" L x 24" W x 53" H / 147cm x 60cm x 134cm
Color	Frame: Charcoal Cover: Cool Gray

WARRANTY

Parts	3 years
Labor	1 year

ENHANCEMENTS

Fit-Key® software	Optional Fit-Key® software makes it easy to create specialized programs and track data
Entertainment options	Yes

ACCESSORIES

Adjustable Pedal Cranks	Provides three adjustments (5", 6" and 7") (12cm, 15cm, 17cm) to give the user versatility to adjust their range of motion
Low Support Boots	Keeps the foot on the pedal with sturdy straps around the back, top, and front of the foot.
Sports Performance Pedals	High-quality design for high-intensity training. Compatible with Shimano SPD cycling shoes. Reversible for cleat connection. Easily removable.
Wheelchair platform	Provides a sturdy foundation for wheelchair use
Heavy duty wheelchair ramps	For manual and powered wheelchair access. Provides increased stabilization
Polar Chest Strap	Transmits heart rate directly to console - works on all SCIFIT products
Seat Belt	Adjustable seat belt for added security and stability. Easy to attach and remove. Storage pouch included. Fits Premium seat only.



SCIFIT Systems, Inc. reserves the right to change or alter specifications at any time without notice and without incurring any obligation.

©Copyright 2015, SCIFIT Systems, Inc. SCIFIT, Fit-Quik, and Fit-Key, are registered trademarks of SCIFIT Systems, Inc. Intelli-Fit and Prescribed for Progress are trademarks of SCIFIT Systems, Inc. Fit-Key Technology is Patented.

ISO7000R Specs - ISO7011R-INT - 2015 - v1