

Identity: Kold Wrap (Cool Bandage) Compression Product - Part #'s 0007, 0248, 0372, 0444, 1632, 0372H

Section I

Manufacturer: Romaine, Inc., d.b.a. Koldcare	Information Telephone: 800-294-7101
Address: 2026 Sterling Ave. Elkhart, Indiana 46516	Date Prepared: 6/28/2002

Section II - Hazardous Ingredients/Identity Information

Hazardous Components: None

Section III - Significant Components

Gel	Carrier Material
Water Polyvinyl Alcohol Sodium Borate (< 1%) Preservatives	Polyurethane Foam

Section IV - Physical/Chemical Characteristics

Boiling Point: 212° F	Viscosity: High
Melting Point: Not Determined	Gel Color: Clear
Flash Point: >200° F	Gel Odor: None

Section V - Fire and Explosion Hazard Data

Flash Point: N/A	Reactivity: None
Extinguishing Media: This is a nonflammable gel solution	Hazardous Waste Defined By 40 CFR 261.4: No

Section VI - Reactivity Data

Stability: Yes	Hazardous Decomposition: None
Incompatibility: None Known	Conditions to Avoid: Temperatures below 32° F and above 90° F

Section VII - Health Hazard Data

Oral Toxicity: Non-toxic by CPSC, OSHA and DOT Toxicity Criteria	Dermal Irritation: Non-irritating according Primary Skin Irritation in Humans
--	---

Emergency First Aid Procedures: Notes to physician, in the unlikely event of rapid ingestion of large volumes, observe patient for diarrhea. (Registered with the Poison Control Center of Indiana - 317-630-6382)

Section VII - Other

No spill procedures, disposal requirements, special protection or precautions, and no DOT shipping regulations.