

Material Safety Data Sheet	<b>North Coast Ultrasound Gel</b>
Manufacturer	North Coast Medical, Inc. 18305 Sutter Blvd. Morgan Hill, CA 95037
Emergency Phone #	(408) 776-5088

HAZARDOUS INGREDIENTS				
Ingredients	%	CAS Number	LD/50 (Species & Route)	LC/80 (Species & Route)
THIS PRODUCT IS KNOWN TO BE A NON-HAZARDOUS PRODUCT				

PHYSICAL DATA			
Physical State	liquid	Vapor Pressure	n/a
Density (g/cm <sup>3</sup> )	1.03	Vapor Density	n/a
Boiling Point (°C)	100	Percent Volatile by Volume	n/a
Freezing Point (°C)	0	Evaporation Rate	n/a
PH	6.50-7.50	Solubility in Water (20°C)	good
Odor & Appearance	blue viscous gel		
Odor Threshold	n/a		

FIRE AND EXPLOSION HAZARD DATA	
Flammability Classification	non-flammable
Flammable Limits	
Lower Explosion Level (% by Volume)	n/a
Upper Explosion Level (% by Volume)	n/a
Flashpoint	n/a
Extinguish With	n/a

REACTIVITY DATA	
Chemical Stability	n/a
Incompatible Materials	n/a
Hazardous Polymerization	n/a

Hazardous Decomposition	n/a
-------------------------	-----

TOXICOLOGICAL PROPERTIES	
ROUTE OF ENTRY	SYMPTOMS
Skin Contact	n/a
Skin Absorption	n/a
Eye Contact	may cause redness to eye
Inhalation Acute	n/a
Inhalation Chronic	n/a
Ingestion	may cause stomach upset if ingested
Effects of Chronic Exposure to Product	n/a
Irritancy of Product	n/a
Exposure Limits of Product	n/a
Sensitization to Product	n/a
Synergistic Materials	n/a

PREVENTATIVE MEASURES	
Engineering Controls	n/a
Waste Disposal	dispose of in accordance to local, provincial and federal regulations.
Leak and Spill Procedure	scoop up excess, clean with appropriate cleaner and rinse with water.
Handling Procedures & Equipment	no special requirements
Storage Requirements	store away from heat and extreme cold
Special Shipping Instructions	should not freeze

FIRST AID MEASURES	
Skin Contact	n/a
Eye Contact	flush with water for 15 minutes while holding eyelids open. Seek medical attention if irritation still exists.
Inhalation	n/a
Ingestion	induce vomiting. Seek medical attention.

PREPARATION INFORMATION	
Date	January 1, 2000