

LITE-TOUCH PEN

Reduces hand and joint fatigue while writing



Easy To Hold, Easy To Write

- Requires just minimal pressure to write.
- Entire pen is built-up to provide a comfortable hold.
- Slight indentation in thumb area aids in gripping.
- Made of lightweight plastic with a clip top.
- Medium point, black ink.
- Use Parker #30316 ink cartridge for replacement.

*Write
comfortably
Requires very
little pressure!*

The Lite-Touch Pen requires almost no pressure to write. The built-up, textured surface makes it easier to hold than standard size pens. A gentle indentation on the shaft helps relieve painful thumbs and aids in gripping. This hollow, lightweight pen is ideal for persons with arthritis.

Circumference measures 2¼" (6cm). Medium point, black ink. Replace ink cartridge with Parker #30316 refill.

NC21020



North Coast Medical, Inc.
8100 Camino Arroyo • Gilroy, CA 95020
Local/Int'l: 408-776-5000 • Toll-Free: 800-821-9319
Toll-Free Fax: 877-213-9300 • www.ncmedical.com



REV21020_0611