

This economical orthosis has been designed for optimal fit and comfort.



### Norco™ D-Ring Wrist Orthosis

Left	Right	
NC15800-0	NC15800-1	X-Small
NC15800-2	NC15800-3	Small
NC15800-4	NC15800-5	Medium
NC15800-6	NC15800-7	Large
NC15800-8	NC15800-9	X-Large

### Norco™ D-Ring Wrist Orthosis (Long)

Left	Right	
NC15829-0	NC15829-1	X-Small
NC15829-2	NC15829-3	Small
NC15829-4	NC15829-5	Medium
NC15829-6	NC15829-7	Large
NC15829-8	NC15829-9	X-Large

### Norco™ D-Ring Black Wrist Orthosis

Left	Right	
NC16150	NC16151	X-Small
NC16152	NC16153	Small
NC16154	NC16155	Medium
NC16156	NC16157	Large
NC16158	NC16159	X-Large

The **Norco™ D-Ring Wrist Orthosis** is used to treat wrist strains, sprains, arthritis, carpal tunnel syndrome and tendinitis. The comfortable, cotton/polyester fabric blend allows full dexterity of the fingers and thumb yet limits wrist motion. Enlarged thumb hole allows full thumb use, even if swollen. The metal volar stay can be removed and adjusted for an optimal fit. The double layer of sewn cotton stockinette helps cushion the strap closure. The D-ring straps allow quick and easy application and tension control. Machine wash and air dry. Latex-free.

The Norco™ Wrist Orthosis measures approximately 6 ½" (16 cm) in length.

The Norco™ Long Wrist Orthosis measures approximately 8 ½" (21 cm) in length.

### Norco™ D-Ring Wrist Orthosis

Measure the circumference at the wrist crease to determine size.

	Wrist Circumference
X-Small	Up to 5¾" (11 to 14cm)
Small	5" to 6" (14 to 17cm)
Medium	6" to 7" (17 to 18cm)
Large	7" to 8" (18 to 20cm)
X-Large	Over 8" (Over 20cm)

# D-Ring Wrist/Thumb Orthoses

Comfortable, padded cotton support for thumb and wrist injuries.



The **Norco™ D-Ring Wrist/Thumb Orthosis** is ideal for cumulative trauma injuries, de Quervain's and tendinitis. The Norco™ D-Ring Wrist/Thumb Orthosis helps decrease pain and inflammation in the wrist and thumb due to arthritis, sprains, strains, joint instabilities, and post-surgical repair. This soft, cotton/polyster blend is circumferentially designed with a double-layer stockinette under the D-ring strap closures for comfort and easy application. It features three straps and a wraparound thumb piece that accommodates all thumb sizes. Includes a metal volar wrist stay and flexible plastic stays along the dorsal opening and along the radial thumb for additional stability. The removable metal stay can be bent to conform as needed. Machine wash and air dry. Latex-free.

The Norco™ Wrist/Thumb Orthosis measures approximately 6½" (16 cm) in length.

## Norco™ D-Ring Wrist Orthosis

Left	Right	
NC15825-01	NC15825-02	<b>X-Small</b>
NC15825-03	NC15825-04	<b>Small</b>
NC15825-05	NC15825-06	<b>Medium</b>
NC15825-07	NC15825-08	<b>Large</b>
NC15825-09	NC15825-10	<b>X-Large</b>

## Norco™ D-Ring Wrist/Thumb Orthosis

Measure the circumference at the wrist crease to determine size.

### Distal Palmar Crease Circumference

<b>X-Small</b>	6" to 7"	(15 to 18cm)
<b>Small</b>	7" to 7¾"	(18 to 20cm)
<b>Medium</b>	7¾" to 8¾"	(20 to 21cm)
<b>Large</b>	8¾" to 9"	(21 to 23cm)
<b>X-Large</b>	Over 9"	(Over 23cm)

