

Freedom Stim™ Dysphagia Electrodes



Cost-effective, clinically-accepted alternative.

Freedom Stim™ Dysphagia Electrodes are FDA cleared for placement on the anterior portion of the neck.

These replacement electrodes are compatible with existing NMES units designed to treat dysphagia (adjustable wave form: 0-25 mA output).

The extremely low impedance allows for increased patient comfort while using a lower mA setting for a safer, more effective treatment. The higher posts allow good adherence to the lead wires during treatment.

Each electrode has two 2 cm. diameter conductor pads that are made with 100% pure silver bonded to a conductive carbon film, which provide uniform current dispersion. The electrode's flexible, medical-grade, tan foam topcoat contours to the neck area and stays securely in place due to the superior adhesion of the electrode's dual-layered hydrogel.

Each metallized, re-sealable package includes: 2 butterfly electrodes (with dual conductor pads), 3 alcohol skin barrier wipes, 3 lotion packets and instructions for use. Electrode measures 2" x 3.5" (5 cm x 9 cm).

Available in boxes of 12, 30, or 50 treatment packets.

Latex free.

This dysphagia treatment should only be done by a certified SLP, OT or PT who has completed the dysphagia training.



Tan Foam	Prof. Price
NC89800-12	159.00
NC89800-30	379.95
NC89800-50	595.00

Scan this QR code
to view a video on
Freedom Stim™
Dysphagia Electrodes



Single Patient, multiple application use.

Warning: Do not apply to broken skin.

REV89800T_0113

Distributed by:
North Coast Medical, Inc.
www.ncmedical.com

