

SCIFIT

AC5000M Treadmill

Scientific Solutions

- Low starting speed at only .1 and ranges up to 12 mph with speed adjustable in only .1 increments.
- Reversible speed goes up to 4 mph and has a -3% decline for downhill training.
- Elevation ranges from -3%-12% adjustable in .5 increments.
- Side handrail switches to adjust speed and elevation allowing the user multiple areas for easy adjustment.
- Maintenance-free, shock absorbing deck. The 2-ply belt with special carbon fibers woven into the fabric and spaced evenly apart to dissipate static electricity.
- 3 HP AC motor with variable frequency torque vector drive system.
- Generous walking surface of 22"x62".
- Dual emergency stop with safety clip and button by the console.
- Extended handrails included.
- 550 lb. user weight capacity.
- Heart Rate compatible with most wireless technology.

Intelli-Fit™ Features and Programs

- 7" LCD display with user friendly interface using tactile keypad, color coded with audible beep.
- Programs: One button Quick Start, Manual, Heart Rate (chest strap must be used), 7 pre-programmed elevation profiles, 7 pre-programmed speed profiles.
- Administrative mode: includes 6 user defined languages, metric or US units of measure and other custom settings.
- Built-in twin Fans: 3 levels of intensity at low, medium and high.
- Automatic cool down after each exercise is completed.
- Download results to USB flash drive. Simply plug your USB key in the Fit-Key slot and your work-out will be copied in a text file for saving and printing.
- Display Feedback: Heart rate, time, speed, elevation, calories, distance in miles, pace, ascent (ft.) and METs.
- CSAFE port and 8-volt DC power port for TV or other controllers.
- USB charger for phones, mp3 players and other electronic devices.

Warranty

3 year parts, 1 year labor with a 5 year frame and drive system warranty. Exceptions: 1 year warranty on: deck, running belt, contact heart rate grips, heart rate receiver/transmitter. Wear items (grips and rubber foot pads) are excluded from warranty.



Extended handrails
included

Fit-Key[®]
Compatible

Technical Specifications

- Overall Dimensions: 89"L x 34" W x 63" H
- Weight of unit: 483 lbs pre packaged
- Shipping Dimensions and weight:
89"x34"x70" 525 lbs (Skid only)
89" X 34" X 70" 525 lbs Boxed
- Maximum user weight: 550 lbs
- 110v standard; 220v optional
- Electrical requirements: dedicated 15 amp circuit, grounded.

Color

Frame: Silver Powder Coated
Covers: Metallic Gray with UV
Protection

Order Code

SF-AC5000M-INT

www.ncmedical.com
800-821-9319