

SCIFIT

SXT7000e² Total Body Elliptical

Scientific Solutions

- Comfortable 17" stride length that aligns the pivot point of the arm handles with the knee to provide spinal alignment. This provides proper posture and correct biomechanical positioning.
- Optimized torso rotation and core recruitment with 22" arm extension and flexion.
- Oversized orthopedic footbeds with safety edge for optimized comfort and stability.
- Patented Bio-Flex technology moves with the feet helps provide better circulation by eliminating hot spots and numbness and helps improve overall posture and balance by keeping the foot moving.
- Bi-directional resistance - the user can exercise in forward and reverse motion.
- Low starting resistance at only 6 watts.
- 1-20 levels adjustable in .1 increments providing 200 levels of resistance.
- Step-up platform and extended handrails for easy entry/exit and greater support
- Heart Rate compatible with most wireless technology.

Intelli-Fit Features and Programs

- 7" LCD display with user friendly interface using tactile keypad, color coded with audible beep.
- Programs: One button Quick Start, Manual, Hill Profiles, Random Profiles, Heart Rate, Iso-Strength, Constant Work, Heart-Fit Test, Heart-Fit Training (4 levels with 3 programs at each level), Power-Fit Test, Power-Fit Training (4 levels with 3 programs at each level), Stress Test and Fit-Quik™.
- Administrative mode: includes 6 user defined languages, metric or US units of measure and other custom settings.
- Built-in twin Fans: 3 levels of intensity at low, medium and high.
- Power Burst - 8 second burst of Iso-Strength. Available on all programs except: Heart-Fit Test/Training, Power-Fit Test/Training, Heart Rate, Iso-Strength and Stress Test.
- Automatic cool down after each exercise is completed.
- Download results to USB flash drive. Simply plug your USB key in the Fit-Key slot and your work-out will be copied in a text file for saving and printing.
- Display Feedback: Heart rate, time, RPM, calories, distance/miles, level, watts and MET's with two change view screens
- CSAFE port and 8-volt DC power port for TV or other controllers.
- USB charger for phones, mp3 players and other electronic devices.



Fit-Key®
Compatible

Technical Specifications

- Overall Dimensions: 75"L x 30" W x 63" H
- Weight of unit: 336 lbs pre packaged
- Shipping Dimensions and weight:
68"x32"x66" 380 lbs (Skid only)
68 5/8" X 32 1/2" X 67 9/16" 415 lbs Boxed
- Maximum user weight: 425 lbs
- Self Generating with auto recharge battery backup.

Color

Frame: Silver Powder Coated
Covers: Metallic Gray with UV Protection

Warranty

3 year parts, 1 year labor. Exceptions: 1 Year warranty on footbeds and Bio-Flex™ pads. Wear items (rubber hand grips) are excluded from warranty.

Order Code

SF-SXT7000E-INT

www.ncmedical.com
800-821-9319