

SCIFIT

**TC1000
Climber**

Scientific Solutions

- Quiet and smooth movement features independent step action.
- Wide variety of user-friendly programs that allow users to change their climbing speed at any time in one-foot-per-minute increments.
- 1-20 levels adjustable in .1 increments providing 200 levels of resistance.
- Self generating with optional AC adapter.
- Large, orthopedic footbeds with safety edge for superior traction and comfort.
- Ergonomic handlebars with contact heart rate permit multiple hand positions.
- Climbing speed 5' to 110'.
- Step range at 10.5".
- Heart Rate compatible with most wireless technology.

Intelli-Fit® Features and Programs

- 7" LCD display with user friendly interface using tactile keypad, color coded with audible beep.
- Programs: One button Quick Start, Manual, 7 Hill Profiles, unlimited number of Random Profiles, Heart Rate, Fit Quik®.
- Administrative mode: includes 6 user defined languages, metric or US units of measure and other custom settings.
- Built-in twin Fans: 3 levels of intensity at low, medium and high.
- Automatic cool down after each exercise is completed.
- Download results to USB flash drive. Simply plug your USB key in the Fit-Key slot and your work-out will be copied in a text file for saving and printing.
- Display Feedback: Heart rate, time, step, steps per minutes (SPM), calories, ascent (ft), and level.
- CSAFE port and 8-volt DC power port for TV or other controllers.
- USB charger for phones, mp3 players and other electronic devices.



Fit-Key®
Compatible

Technical Specifications

- Overall Dimensions: 39" L x 26" W x 70" H
- Weight of unit: 182 lbs pre packaged
- Shipping Dimensions and weight:
68" x 32" x 66" 227 lbs (Skid only)
68 5/8" X 32 1/2" X 67 9/16" 262 lbs Boxed
- Maximum user weight: 375 lbs
- Minimum user weight: 100 lbs
- Self Generating with auto recharge battery backup. Low voltage AC adapter included.

Color

Frame: Silver Powder Coated
Covers: Metallic Gray with UV Protection

Warranty

3 year parts, 1 year labor. Exceptions: 1 Year warranty on footbeds. Wear items (rubber hand grips) are excluded from warranty.

Order Code

SF-TC1002-INT
www.ncmedical.com
800-821-9319